

Product datasheet for TP308985

WDR92 (NM_138458) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human WD repeat domain 92 (WDR92), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208985 protein sequence Red =Cloning site Green =Tags(s)
	MSAFEKPKQIIAHIQKGFNYTVFDCKWVPCSAKFVTMGNFARGTGVIQLYEIQHGDLKLLREIEKAKPIKC GTFGATSLQQRYLATGDFGGNLHIWNLEAPEMPVYSVKGHKEIINAIDGIGGLGIGEGAPEIVTGSRDGT VKVWDPRQKDDPVANMEPVQGENKRDCWTVAFGNAYNQEERVVCAGYDNGDIKFLDLRNMALRWETNIKN GVCSLEFDRKDISMNKLVATSLEGKFHVFD MRTQHPTKGFASVSEKAHKSTVWQVRHLPQNRELFLTAGG AGGLHLWKYEYPIQRSKKDSEGIEMGVAGSVSLLQNVTLSTQPISSLDWSPDKRGLCVCSFDQTVRVLI VTKLNKI
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	39.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_612467
Locus ID:	116143



[View online »](#)

UniProt ID: [Q96MX6](#), [A0A140VK67](#)

RefSeq Size: 2280


Cytogenetics: 2p14

RefSeq ORF: 1071

Summary: This gene encodes a protein with two WD40 repeat domains thought to be involved in an apoptosis via activation of caspase-3. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2012]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified WDR92 protein (Cat# TP308985). The protein was produced from HEK293T cells transfected with WDR92 cDNA clone (Cat# [RC208985]) using MegaTran 2.0 (Cat# [TT210002]).