

Product datasheet for **TP308988M**

EEPD1 (NM_030636) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human endonuclease/exonuclease/phosphatase family domain containing 1 (EEPD1), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC208988 protein sequence Red =Cloning site Green =Tags(s)
	MGSTLGCHRSIPRDPDDLHSHSRKFSAACNFSNVLNQRNLNINTATEEELMTLPGVTRAVARSIVEYREY IGGFKKVEDLALVSGVGATKLEQVKFEICVSSKGSSAQHSPSSLRDLLAEQQPHHLATAVPLTPRVNIN TATPAQLMSVRGLSEKMALSIQDFRREHGPFRSVEDLVMDGINAAFLDRIRHQVFAERSRPPSTHTNGG LTFTAKPHPSPTSLSLQSEDLPLPPGGPTQIISTRPSVEAFGGTRDGRPVLRATWNLQGCSVEKANNPG VREVCMITLLENSIKLLAVQELLDRLEALEKFCTELNQPTLPNIRKWKGPGRGCWKAVVAEKPSNQLQKGAG YAGFLWDAAAGMELRDAGSQESSPSNGHGKLAGPSPYLGRFKVGSHDLTLVNLHLAALTLLGSENPSKN H SDGHRLASFAQTLQETLKGEKDVILGDFGQGPDSNDYDILRKEKFHHLIPAHTFTNISTKNPQGSKSLD NIWISKSLKKVFTGHWAVVREGLTNPWIPDNWSWGGVASEHCPVLAEFYTEKDWSKKDAPRNGSGVALE R SEANIKHER SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	62.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_085139
Locus ID:	80820
UniProt ID:	Q7L9B9
RefSeq Size:	4687
Cytogenetics:	7p14.2
RefSeq ORF:	1707
Synonyms:	HSPC107

Product images:



Coomassie blue staining of purified EEPD1 protein (Cat# [TP308988]). The protein was produced from HEK293T cells transfected with EEPD1 cDNA clone (Cat# [RC208988]) using MegaTran 2.0 (Cat# [TT210002]).