

Product datasheet for **TP309049M**

THNSL2 (NM_018271) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human threonine synthase-like 2 (<i>S. cerevisiae</i>) (THNSL2), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209049 representing NM_018271 Red =Cloning site Green =Tags(s)
	 MWYVSTRGVAPRVNFEGALFSGYAPDGGLFMPEELPQLDRGTLCQWSTLSYPGLVKELCALFIGSELLPK DELNDLIDRAFSRFRHREVHLSRLRNLNVLELWHGVTYAFKDLSSLCTTQFLQYFLEKREKHVTWVG TSGDTGSAAIESVQGAKNMDIIVLLPKGHCTKIQELQMTTVLKQNVHVFVGVGNSDELDEPIKTVFADVA FVKKHNLMSLNSINWSRVLVQMAHHFFAYFQCTPSLDTHPLPLVEVVVPTGAAGNLAAGYIAQKIGLPIR LVVAVNRNDIIHRTVQQGDFSLSEAVKSTLASAMDIQVPYNMERVFWLLSGSDSQVTRALMEQFERTQSV NLPKELHSKLSEAVTSVSVSDEAITQTMGRCWDENQYLLCPHSAVAVNYHYQQIDRQQPSTPRCCLAPAS AAKFPEAVLAAGLTPETPAEIVALEHKETRCTLMRRGDNWMLMLRDTIEDLSRQWRSHALNTSQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	53.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_060741




[View online »](#)

Locus ID: 55258
UniProt ID: [Q86YJ6](#)
RefSeq Size: 2150
Cytogenetics: 2p11.2
RefSeq ORF: 1452
Synonyms: SOFAT; THS2; TSH2

Summary: This gene encodes a threonine synthase-like protein. A similar enzyme in mouse can catalyze the degradation of O-phospho-homoserine to a-ketobutyrate, phosphate, and ammonia. This protein also has phospho-lyase activity on both gamma and beta phosphorylated substrates. In mouse an alternatively spliced form of this protein has been shown to act as a cytokine and can induce the production of the inflammatory cytokine IL6 in osteoblasts. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2011]

Product images:



Coomassie blue staining of purified THNSL2 protein (Cat# [TP309049]). The protein was produced from HEK293T cells transfected with THNSL2 cDNA clone (Cat# [RC209049]) using MegaTran 2.0 (Cat# [TT210002]).