

Product datasheet for **TP309050L**

CCDC25 (NM_018246) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coiled-coil domain containing 25 (CCDC25), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209050 representing NM_018246 Red =Cloning site Green =Tags(s)

MVFYFTSSSVNSSAYTIYMGKDKYENEDLIKHGWPEDIWFHVDKLSSAHVYLRHLKGENIEDIPKEVLMD
CAHLVKANSIQGCKMNNVNVVYTPWSNLKKTADMDVGQIGFHRQKDVKIVTVEKKVNEILNRLEKTKVER
FPDLAAEKECRDREERNEKKAQIQEMKKREKEEMKKKREMDLRSYSSLMKVENMSSNQDGNDSDEFM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	24.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_060716
Locus ID:	55246
UniProt ID:	Q86WR0
RefSeq Size:	3653




[View online »](#)

Cytogenetics: 8p21.1

RefSeq ORF: 624

Summary: Transmembrane receptor that senses neutrophil extracellular traps (NETs) and triggers the ILK-PARVB pathway to enhance cell motility (PubMed:32528174). NETs are mainly composed of DNA fibers and are released by neutrophils to bind pathogens during inflammation (PubMed:32528174). Formation of NETs is also associated with cancer metastasis, NET-DNA acting as a chemotactic factor to attract cancer cells (PubMed:32528174). Specifically binds NETs on its extracellular region, in particular the 8-OHdG-enriched DNA present in NETs, and recruits ILK, initiating the ILK-PARVB cascade to induce cytoskeleton rearrangement and directional migration of cells (PubMed:32528174). In the context of cancer, promotes cancer metastasis by sensing NETs and promoting migration of tumor cells (PubMed:32528174). [UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified CCDC25 protein (Cat# [TP309050]). The protein was produced from HEK293T cells transfected with CCDC25 cDNA clone (Cat# [RC209050]) using MegaTran 2.0 (Cat# [TT210002]).