

Product datasheet for **TP309088L**

DPH4 (DNAJC24) (NM_181706) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human DnaJ (Hsp40) homolog, subfamily C, member 24 (DNAJC24), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209088 protein sequence Red =Cloning site Green =Tags(s) MAVEQMPKKDWYSILGADPSANISDLKQKYQKLILMYHPDKQSTDVPAGTVEECVQKFIEIDQAWKILGN EETKREYDLQRCEDDLRNVGPDVAQVYLEEMSWNEGDHSFYLSRCGGKYSVSKDEAEVSLISCDTCSL IIELLHYN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	17 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_859057</u>
Locus ID:	120526
UniProt ID:	<u>Q6P3W2</u>


[View online »](#)

RefSeq Size:	3000
Cytogenetics:	11p13
RefSeq ORF:	444
Synonyms:	DPH4; JJJ3; ZCSL3
Summary:	Diphthamide is a unique posttranslationally modified histidine found only in translation elongation factor-2 (EEF2; MIM 130610). This modification is conserved from archaeobacteria to humans and serves as the target for ADP-ribosylation and inactivation of EEF2 by diphtheria toxin (DT) and Pseudomonas exotoxin A. DPH4 is 1 of several enzymes involved in synthesis of diphthamide in EEF2 (Liu et al., 2004 [PubMed 15485916]).[supplied by OMIM, Mar 2008]

Product images:



Coomassie blue staining of purified DNAJC24 protein (Cat# [TP309088]). The protein was produced from HEK293T cells transfected with DNAJC24 cDNA clone (Cat# [RC209088]) using MegaTran 2.0 (Cat# [TT210002]).