

Product datasheet for **TP309186**

TTC23L (NM_144725) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human tetratricopeptide repeat domain 23-like (TTC23L), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209186 protein sequence Red =Cloning site Green =Tags(s)
	<p>MQASPIRIPTVSNIDIDWDFCFRMSQQTEIPAHQQTDELYPTGGCGESEETKAKEKEKAIDCMSHPKEKL AQSQKKVAQLIKEKMNTQANKELIRCIVLSRIIFGDHHWKCARALANLAYGYLTRGLPVQAKKHATSAK NTLLTWKANTTSNKEKEEILEALVKLYYTLGVAWLLQNRGREAYFNLQKAERNMKELKELYKGGVCELQV SENDLTALGRASLAHRLNLALAYFEKAIGDVIAAKGDRTSDLISLYEEAAQIEQLRRNHNQAIQYLQQ AHSVCVSLFTEVSPKTAEMSALLAKAYAMSGEAQHRDAVEIYFIRSINAYRATLGSEDFETLSTTEEFCK WLVQNGEKQDK</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	40.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_653326</u>
Locus ID:	153657



[View online »](#)

UniProt ID: [Q6PF05](#)
 RefSeq Size: 1369
 Cytogenetics: 5p13.2
 RefSeq ORF: 1083
 Synonyms: MC25-1

Product images:



Coomassie blue staining of purified TTC23L protein (Cat# TP309186). The protein was produced from HEK293T cells transfected with TTC23L cDNA clone (Cat# [RC209186]) using MegaTran 2.0 (Cat# [TT210002]).