

Product datasheet for **TP309206L**

WNT5A (NM_003392) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human wingless-type MMTV integration site family, member 5A (WNT5A), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209206 protein sequence Red =Cloning site Green =Tags(s)

MKKSIGILSPGVALGMAGSAMSSKFFLVALAIFFSFAQVVEANSWWSLGMNPNVQMSEVYIIGAQPLCS
QLAGLSQGQKKLCHLYQDHMQYIGEGAKTGIKECQYQFRHRRWNCSTVDNTSVFGRVMQIGSRETAFTYA
VSAAGVWNAMSRAAREGELSTCGCSRAARPKDLPRDWLWGGCGDNIDYGYRFAKEFVDARERERIHAKGS
YESARILMNLHNNEAGRRTVYNLADVACKCHGVSGSCSLKTCWLQLADFRKVGDLALKEKYDSAAAMRLNS
RGKLVQVNSRFNSPTTQDLVYIDPSPDYCVRNESGSLGTQGRLCNKTSEGMDGCELMCCGRGYDQFKTV
QTERCHCKFWCCYVKCKKCTEIVDQFVCK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	42.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_003383</u>



[View online »](#)

Locus ID: 7474

UniProt ID: [P41221](#), [A0A384N611](#), [B3KQX9](#)

RefSeq Size: 6194

Cytogenetics: 3p14.3

RefSeq ORF: 1140


Synonyms: hWNT5A

Summary: The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene encodes a member of the WNT family that signals through both the canonical and non-canonical WNT pathways. This protein is a ligand for the seven transmembrane receptor frizzled-5 and the tyrosine kinase orphan receptor 2. This protein plays an essential role in regulating developmental pathways during embryogenesis. This protein may also play a role in oncogenesis. Mutations in this gene are the cause of autosomal dominant Robinow syndrome. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Jan 2012]

Protein Families: Adult stem cells, Cancer stem cells, Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein, Stem cell relevant signaling - Wnt Signaling pathway

Protein Pathways: Basal cell carcinoma, Hedgehog signaling pathway, Melanogenesis, Pathways in cancer, Wnt signaling pathway

Product images:



Coomassie blue staining of purified WNT5A protein (Cat# [TP309206]). The protein was produced from HEK293T cells transfected with WNT5A cDNA clone (Cat# [RC209206]) using MegaTran 2.0 (Cat# [TT210002]).