

Product datasheet for **TP309218L**

Desmocollin 2 (DSC2) (NM_024422) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens desmocollin 2 (DSC2), transcript variant Dsc2a, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209218 representing NM_024422 Red =Cloning site Green =Tags(s)

MEAARPSGSWNGALCRLLLTLAILIFASDACKNVTLHVPSKLD AEKLVGRVNLKECF TAANLIHSSDPD
FQILEDGSVYTTNTILLSSEKRSFTILLSNTENQEKKKIFVLEHQTKVLKKRHTKEKVLRRAKRRWAPI
PCSMLENSLGPFLFLQQVQSDTAQNYTIYYSIRGPGVDQEPNLFYVERDTGNLYCTRPVDREQYESFE
IIAFATTPDGYTPELPLPLIIKIEDENDNYPIFTEETYFTIFENCRVGTTVGQVCATDKDEPDTMHTRL
KYSIIGQVPPSPTLFSMHPTTGVITTTSSQLDRELIDKYQLKIKVQMDMGQYFGLQTTSTCIINIDDVND
HLPTFTRTSYVTSVEENTVDVEILRVTVEDKDLVNTANWRANYTILKGNENGNFKIVTDAKTNEGVL CV
KPLNYEEKQMILQIGVVNEAPFSREASPRSAMSTATVTNVNVEDQDEGPCNPPIQTVRMKENA EVGTT S
NGYKAYDPETRSSSGIRYKKLTDPTGWVTIDENTGSIKVFRLDREAETIKNGIYNITVLASDQGGRTCT
GTLGIIQDVNDNSPFIPKKTVIICKPTMSSAEIVAVDPDEPIHGPPDFSLSTSEVQRMWRLKAIN D
TAARLSYQNDPPFGSYVVPITVRDLGMSSVTSLDVTLDCITENDCTHRVDPRI GGGGVQLGKWAILAI
LLGIALLFILFTLVCGASGTSKQPKVIPDDLAQQNLIVSNTEAPGDDKVYSANGFTTQTVGASAQGVCG
TVGSGIKNGGQETIEMVKGGHQTSESCRAGHHHTLDSCRGGHTEVDNCRYTYSEWHSFTQ PRLGEKVYL
CNQDENHKHAQDYVLTYNYEGRGSVAGSVGCCSERQEEDGLEFLDNLEPKFRTLAEACMKR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	97 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_077740
Locus ID:	1824
UniProt ID:	Q02487
RefSeq Size:	5140
Cytogenetics:	18q12.1
RefSeq ORF:	2703
Synonyms:	ARVD11; CDHF2; DG2; DGII/III; DSC3
Summary:	This gene encodes a member of the desmocollin protein subfamily. Desmocollins, along with desmogleins, are cadherin-like transmembrane glycoproteins that are major components of the desmosome. Desmosomes are cell-cell junctions that help resist shearing forces and are found in high concentrations in cells subject to mechanical stress. This gene is found in a cluster with other desmocollin family members on chromosome 18. Mutations in this gene are associated with arrhythmogenic right ventricular dysplasia-11, and reduced protein expression has been described in several types of cancer. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2015]
Protein Families:	Transmembrane
Protein Pathways:	Arrhythmogenic right ventricular cardiomyopathy (ARVC)

Product images:



Coomassie blue staining of purified DSC2 protein (Cat# [TP309218]). The protein was produced from HEK293T cells transfected with DSC2 cDNA clone (Cat# [RC209218]) using MegaTran 2.0 (Cat# [TT210002]).