

Product datasheet for TP309226L


ADAMTS4 (NM_005099) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ADAM metallopeptidase with thrombospondin type 1 motif, 4 (ADAMTS4), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209226 representing NM_005099 Red =Cloning site Green =Tags(s) MSQTGSHPGRGLAGRWLWGAQPCLLLPIVPLSWLVWLLLLLLASLLPSARLASPLPREEEIVFPEKLNGS VLPGSGAPARLLCRLQAFGETLLLELEQDSGVQVEGLTVQYLGQAPELLGGAEPGYLTGTINGDPESVA SLHWDGGALLGVLQYRGAEHLQPLEGGTPNSAGGPGAHLRRKSPASGQGPMCNVKAPLGSPSPRRR A KRFASLSRFVETLVVADDKMAAFHGAGLKRYLLTVMAAAAKAFKHPSIRNPVSLVTRLVILGSGEEGPQ VGPSAAQTLRSFCWQRGLNTPENSDDPHFDTAILFTRQDLCGVSTCDTLGMADVGTVCDPARSCAIVED DGLQSAFTAAHELGHVFNVLHDNSKPCISLNGPLSTSRHVMAPVMAHVDPEEPWSPCSARFITDFLDNG Y GHCLLDKPEAPLHLPVTFPGKDYDADRQCQLTFGPDSRHCPQLPPPCAALWCSGHLNGHAMCQTKHSP WA DGTPCGPAQACMGGRCLHMDQLQDFNIPQAGGWGPWGPWGDCSRTC GGGVQFSSRDCTRPVTRNG GKYCE GRRARFRSCNTEDCPTGSALTFREEQCAAYNHRTDLFKSFPGMDWVPRYTGAVPDQCKLTCQAQALG Y YYVLEPRVVDGTPCSPDSSSVQVQGRCIHAGCDRIIGSKKKFDKCMVCGGDGSGCSKQSGSFRKFRYGYN NVVTIPAGATHILVRQQGNPGHRSIYLALKLPDGSYALNGEYTLMPSPTDVVLPGAVSLRYSGATAASET LSGHGPLAQPLTLQVLVAGNPQDTRLRYSFFVPRPTPSTPRPTPQDWLHRRRAQILEILRRRPWAGRK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	90 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining


[View online »](#)

Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_005090</u>
Locus ID:	9507
UniProt ID:	<u>O75173</u>
RefSeq Size:	4342
Cytogenetics:	1q23.3
RefSeq ORF:	2511
Synonyms:	ADAMTS-2; ADAMTS-4; ADMP-1
Summary:	<p>This gene encodes a member of the ADAMTS (a disintegrin and metalloproteinase with thrombospondin motifs) protein family. Members of this family share several distinct protein modules, including a propeptide region, a metalloproteinase domain, a disintegrin-like domain, and a thrombospondin type 1 (TS) motif. Individual members of this family differ in the number of C-terminal TS motifs, and some have unique C-terminal domains. The enzyme encoded by this gene lacks a C-terminal TS motif. The encoded preproprotein is proteolytically processed to generate the mature protease. This protease is responsible for the degradation of aggrecan, a major proteoglycan of cartilage, and brevican, a brain-specific extracellular matrix protein. The expression of this gene is upregulated in arthritic disease and this may contribute to disease progression through the degradation of aggrecan. Alternative splicing results in multiple transcript variants, at least one of which encodes an isoform that is proteolytically processed. [provided by RefSeq, Feb 2016]</p>
Protein Families:	Druggable Genome, Protease, Secreted Protein, Transmembrane

Product images:


Coomassie blue staining of purified ADAMTS4 protein (Cat# [TP309226]). The protein was produced from HEK293T cells transfected with ADAMTS4 cDNA clone (Cat# [RC209226]) using MegaTran 2.0 (Cat# [TT210002]).