

Product datasheet for TP309273L

Cullin 5 (CUL5) (NM_003478) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens cullin 5 (CUL5), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209273 protein sequence Red=Cloning site Green=Tags(s)

MATSNLLKNKGS LQFEDK WDFMRPIVLKLLRQESVTKQQWFDLFSVDHAVCLWDDKGP AKIHQALKEDIL
EFIKQAQARVLSHQDDTALLKAYIVEWRKFFTC DILPKPFCQLEITLMGKQGSNKKSNVEDSIVRKLML
DTWNESIFSNIKNRLQDSAMKLVHAERLGEAFDSQLVIGVRESYVNLCSNPEDKLQIYRDNFEKAYLDST
ERFYRTQAPSYLQQNGVQNYMKYADAKLKEEEKRALRYLETRRECNSVEALMECCVNALVTSFKETILAE
CQGMIKRNETEKLHLMFSLMDKVPNGIEPMLKDLEEHIISAGLADMVAAAETITTDSEKYVEQLLTLFNR
FSKLVKEAFQDDPRFLTARDKAYKAVVNDATIFKLELPLKQKGVGLKTQPESKCEPELLANYCDMLLRKTP
LSKKL TSEEIEAKLKEVLLVLKYVQNKDVMRYHKAHLTRRLIDISADSEIEENMV EWLREVGMPADYV
NKLARMFQDIKVS EDLNQAFKEMHKNNK LALPADSVNIKLNAGAWSRSSEKVFVSLPTELEDL IPEVEE
FYKKNHSGRKLHW HHLMSNGIITFKNEVGQYDLEVTTFQLAVLFAWNQRPREKISFENLKLATELPDAEL
RRTLWLSLVAFPKLRQVLLYEPQVNSPKDFTEGTLFSVNQEFSLIKNAKVQKRGKINLIGRLQLTTERMR
EEENEGIVQLRILRTQEAI IQIMKMRKKISNAQLQTELVEILKNMFLPQKKMIKEQIEWLIEHKYIRRDE
SDINTFIYMA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	90.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_003469
Locus ID:	8065
UniProt ID:	Q93034
RefSeq Size:	6408
Cytogenetics:	11q22.3
RefSeq ORF:	2340
Synonyms:	CUL-5; VACM-1; VACM1
Summary:	Core component of multiple SCF-like ECS (Elongin-Cullin 2/5-SOCS-box protein) E3 ubiquitin-protein ligase complexes, which mediate the ubiquitination and subsequent proteasomal degradation of target proteins. As a scaffold protein may contribute to catalysis through positioning of the substrate and the ubiquitin-conjugating enzyme. The functional specificity of the E3 ubiquitin-protein ligase complex depends on the variable substrate recognition component. ECS(SOCS1) seems to direct ubiquitination of JAK2. Seems to be involved in proteasomal degradation of p53/TP53 stimulated by adenovirus E1B-55 kDa protein. May form a cell surface vasopressin receptor.[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome
Protein Pathways:	Ubiquitin mediated proteolysis

Product images:



Coomassie blue staining of purified CUL5 protein (Cat# [TP309273]). The protein was produced from HEK293T cells transfected with CUL5 cDNA clone (Cat# [RC209273]) using MegaTran 2.0 (Cat# [TT210002]).