

Product datasheet for **TP309280**

TMEM88 (NM_203411) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human transmembrane protein 88 (TMEM88), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209280 representing NM_203411 Red =Cloning site Green =Tags(s) MADVPGAQRAVPGDGPEPRDPLDCWACAVLTAQNLLVAAFNLLLLVLVLTILLPAVTMLGFGFLCHSQ FLRSQAPPCTAHLRDPGFTALLVTGFLLLVPLLVLALASYRRRLCLRLRLADCLVPYSRALYRRRRAPQPR QIRASPGSQAVPTSGKVVV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	17.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_981956
Locus ID:	92162
UniProt ID:	Q6PEY1
RefSeq Size:	887



[View online »](#)


Cytogenetics: 17p13.1

RefSeq ORF: 477

Summary: Inhibits the Wnt/beta-catenin signaling pathway. Crucial for heart development and acts downstream of GATA factors in the pre-cardiac mesoderm to specify lineage commitment of cardiomyocyte development.[UniProtKB/Swiss-Prot Function]

Protein Families: Transmembrane

Product images:



Coomassie blue staining of purified TMEM88 protein (Cat# TP309280). The protein was produced from HEK293T cells transfected with TMEM88 cDNA clone (Cat# [RC209280]) using MegaTran 2.0 (Cat# [TT210002]).