

Product datasheet for TP309305

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

SOCS3 (NM_003955) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human suppressor of cytokine signaling 3 (SOCS3), 20 μg

Species: Human
Expression Host: HEK293T

Expression cDNA Clone >RC209305 protein sequence or AA Sequence: Red=Cloning site Green=Tags(s)

MVTHSKFPAAGMSRPLDTSLRLKTFSSKSEYQLVVNAVRKLQESGFYWSAVTGGEANLLLSAEPAGTFLI RDSSDQRHFFTLSVKTQSGTKNLRIQCEGGSFSLQSDPRSTQPVPRFDCVLKLVHHYMPPPGAPSFPSPP TEPSSEVPEQPSAQPLPGSPPRRAYYIYSGGEKIPLVLSRPLSSNVATLQHLCRKTVNGHLDSYEKVTQL

PGPIREFLDQYDAPL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 24.6 kDa

Concentration: $>0.05 \mu g/\mu L$ as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by

conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 003946

Locus ID: 9021

UniProt ID: O14543, Q6Fl39





RefSeq Size: 2737

Cytogenetics: 17q25.3

RefSeq ORF: 675

Synonyms: ATOD4; CIS3; Cish3; SOCS-3; SSI-3; SSI3

Summary: This gene encodes a member of the STAT-induced STAT inhibitor (SSI), also known as

suppressor of cytokine signaling (SOCS), family. SSI family members are cytokine-inducible negative regulators of cytokine signaling. The expression of this gene is induced by various cytokines, including IL6, IL10, and interferon (IFN)-gamma. The protein encoded by this gene

can bind to JAK2 kinase, and inhibit the activity of JAK2 kinase. Studies of the mouse counterpart of this gene suggested the roles of this gene in the negative regulation of fetal


liver hematopoiesis, and placental development. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Protein Pathways: Adipocytokine signaling pathway, Insulin signaling pathway, Jak-STAT signaling pathway, Type

II diabetes mellitus, Ubiquitin mediated proteolysis

Product images:



Coomassie blue staining of purified SOCS3 protein (Cat# TP309305). The protein was produced from HEK293T cells transfected with SOCS3 cDNA clone (Cat# [RC209305]) using

MegaTran 2.0 (Cat# [TT210002]).