

Product datasheet for **TP309358M**

GRP (LGALSL) (NM_014181) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human galectin-related protein (HSPC159), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209358 protein sequence Red =Cloning site Green =Tags(s)
	MAGSVADSDAVVKLDDGHLNNSLSSPVQADVYPRLIVPFCGHIKGGMRPGKKVLMGIVDLNPESFAIS LTCGDESDPPADVAIELKAVFTDRQLLRNSCISGERGEEQSAIPYFPFIPDQPFPRVEILCEHPRFRVFVD GHQLFDYHRIQTLSAIDTIKINGDLQITKLG
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	18.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_054900
Locus ID:	29094
UniProt ID:	Q3ZCW2
RefSeq Size:	3842



[View online »](#)

Cytogenetics:	2p14
RefSeq ORF:	516
Synonyms:	GRP; HSPC159
Summary:	Does not bind lactose, and may not bind carbohydrates.[UniProtKB/Swiss-Prot Function]

Product images:

Coomassie blue staining of purified LGALS1 protein (Cat# [TP309358]). The protein was produced from HEK293T cells transfected with LGALS1 cDNA clone (Cat# [RC209358]) using MegaTran 2.0 (Cat# [TT210002]).