

Product datasheet for TP309371

AADACL2 (NM_207365) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human arylacetamide deacetylase-like 2 (AADACL2), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209371 protein sequence Red =Cloning site Green =Tags(s)

MPDNIEESWKIMALDAIAKTCTFTAMCFENMRIMRYEEFISMIFRLDYTQPLSDEYITVTDTTFVDIPVR
LYLPKRKSETRRRAVIYFHGGGFCFGSSKQRAFDLNRWTANTLDAVWVGVDYRLAPQHHFPAQFEDGLA
AVKFFLLEKILTKYGVDPTRICISGDSSGGNLATAVTQQVQND AEIKHKIKMQVLLYPGLQITDSYLP SH
RENEHGIVLTRDVAIKLVSLYFTKDEALPWAMRRNQHMPLESRHLFKFVNWSILLPEKYRKDYVYTEPIL
GGLSYSLPGLTDSRALPLLANDSQLQNLPLTYILTCQHDLRLDDGLMYVTRLRNVGVQVWHEHIEDGIHG
ALSFMTSPFYLRRLGLRIRDMYVSWL DKNL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	43.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_997248
Locus ID:	344752



[View online »](#)

UniProt ID: [Q6P093](#)
RefSeq Size: 1513
Cytogenetics: 3q25.1
RefSeq ORF: 1137

Product images:



Coomassie blue staining of purified AADACL2 protein (Cat# TP309371). The protein was produced from HEK293T cells transfected with AADACL2 cDNA clone (Cat# [RC209371]) using MegaTran 2.0 (Cat# [TT210002]).