

## Product datasheet for **TP309433**

### LZIC (NM\_032368) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human leucine zipper and CTNNBIP1 domain containing (LZIC), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209433 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MASRGKTETSKLKQNLEEQLDRLMQQLQDLEECREELDTDEYEETKKETLEQLSEFNDSLKKIMSGNMTL VDELSGMQLAIQAAISQAFKTPREVIRLFAKKQPGQLRTRLAEMDRDLMVGKLERDLYTQQKVEILTALRK LGEKLTADDEAFLSANAGAILSQFEKVSTDLGSGDKILALASFEVEKTKK
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	21.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_115744</a></u>
Locus ID:	84328
UniProt ID:	<u><a href="#">Q8WZA0</a></u> , <u><a href="#">A0A024R4I7</a></u>
RefSeq Size:	1484




[View online »](#)

Cytogenetics: 1p36.22

RefSeq ORF: 570

### Product images:



Coomassie blue staining of purified LZIC protein (Cat# TP309433). The protein was produced from HEK293T cells transfected with LZIC cDNA clone (Cat# [RC209433]) using MegaTran 2.0 (Cat# [TT210002]).