

## Product datasheet for TP309490M


### C1orf198 (NM\_032800) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 1 open reading frame 198 (C1orf198), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209490 representing NM_032800 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)  MASMAAAIAASRSVMSGNRPLDDRERKRFTYFSSLSPMARKIMQDKEKIREKYGPEWARLPPAQQDEII DRCLVGPRAPAPRDPGDSEELTRFPLRGPTGQKVVRFGDEDLTWQDEHSAPFSWETKSQMEFSISALSI QEPSNGTAASEPRPLSKASQGSQALKSSQGSRSSLDALGPTRKEEEASFWKINAERSRGEGPEAEFQSL TPSQIKSMEKGEKVLPPCYRQEPAPKDREAKVERPSTLRQEQRPLPNVSTERERPQPQAFSSALHEAAP SQLEGKLPSPDVRQDDGEDTLFSEPKFAQVSSSNVWLKTGFDFLDNW  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	36.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_116189</a>
Locus ID:	84886


[View online »](#)

UniProt ID: Q9H425  
 RefSeq Size: 3733  
 Cytogenetics: 1q42.2  
 RefSeq ORF: 981

**Product images:**


Coomassie blue staining of purified C1orf198 protein (Cat# [TP309490]). The protein was produced from HEK293T cells transfected with C1orf198 cDNA clone (Cat# [RC209490]) using MegaTran 2.0 (Cat# [TT210002]).