

Product datasheet for TP309502

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

ASB7 (NM_024708) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human ankyrin repeat and SOCS box-containing 7 (ASB7), transcript

variant 1, 20 µg

Species: Human Expression Host: HEK293T

Expression cDNA >RC209502 protein sequence **Clone or AA Sequence:** Red=Cloning site Green=Tags(s)

MLHHHCRRNPELQEELQIQAAVAAGDVHTVRKMLEQGYSPNGRDANGWTLLHFSAARGKERCVRVFLEHG ADPTVKDLIGGFTALHYAAMHGRARIARLMLESEYRSDIINAKSNDGWTPLHVAAHYGRDSFVRLLLEFK AEVDPLSDKGTTPLQLAIIRERSSCVKILLDHNANIDIQNGFLLRYAVIKSNHSYCRMFLQRGADTNLGR LEDGQTPLHLSALRDDVLCARMLYNYGADTNTRNYEGQTPLAVSISISGSSRPCLDFLQEVTSM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 30.5 kDa

Concentration: >0.05 μg/μL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional

chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling

conditions. Avoid repeated freeze-thaw cycles.

RefSeq: <u>NP 078984</u> **Locus ID:** 140460





RefSeq ORF:

UniProt ID: Q9H672

1769 RefSeq Size:

Cytogenetics: 15q26.3

Summary: The protein encoded by this gene belongs to a family of ankyrin repeat proteins that, along with

four other protein families, contains a C-terminal SOCS box motif. Growing evidence suggests that the SOCS box acts as a bridge between specific substrate-binding domains and the more generic proteins that comprise a large family of E3 ubiquitin protein ligases. In this way, SOCS


box containing proteins may regulate protein turnover by targeting proteins for

polyubiquination and, therefore, for proteasome-mediated degradation. Two alternative transcripts encoding different isoforms have been described. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

822

Product images:



Coomassie blue staining of purified ASB7 protein (Cat# TP309502). The protein was produced from HEK293T cells transfected with ASB7 cDNA clone (Cat# [RC209502]) using MegaTran 2.0 (Cat# [TT210002]).