

Product datasheet for **TP309530L**

DNA Ligase III (LIG3) (NM_013975) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ligase III, DNA, ATP-dependent (LIG3), nuclear gene encoding mitochondrial protein, transcript variant alpha, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209530 protein sequence Red =Cloning site Green =Tags(s)

MSLAFKIFFPQTLRALSRLKELCLFRKHHWRDVRQFSQWSETDLLHGHPLFLRRKPVLSFQGSHLRSRATY
LVFLPGLHVGLCSGPCEMAEQRFCVDYAKRGTAGCKKCKEKIVKGVCRIGKVVNPNPFSSEGGDMKEWYHI
KCMFEKLERARATTKKIEDLTELEGWEELEDNEKEQITQHIADLSSKAAGTPKKKAVVQAKLTTTGQVTS
PVKGASFVTSTNPRKFSGFSAKPNNSGEAPSSPTPKRSLSSSKCDPRHKDCLLREFRKLCAVADNPSYN
TKTQIIQDFLRKGSAGDGFHGDVYLVTKLLPGVIKTVYNLNDKQIVKLFSTRIFNCNPDDMARDLEQGDV
SETIRVFFEQSKSFPPAAKSLLLTIQEVDEFLLRLSKLTKEDEQQALQDIASRCTANDLKCIIRLIKHDH
KMNSGAKHVLDALDPNAYEAFKASRNLDVVERVLHNAQEVEKEPGQRRALSQASLMTPVQPMLEACK
SVEYAMKKCPNGMFSEIKYDGERVQVHKNGDHFSYFSRSLKPVLPKHVAHFKDYIPQAFPGGHSMLDSE
VLLIDNKTGKPLPFGTLGVHKKAAFQDANVCLFVFDYFNDVSLMDRPLCERRKFLHDNMVEIPNRIMF
SEMKRVTKALDLADMTRVIQEGLEGLVLDKVKGTYPGKRHWLKVKKDYLNEGAMADTADLVVLGAFYG
QGSKGGMMSIFLMGCYDPGSQKWCTVTKCAGGHDDATLARLQNELDMVKISKDPSKIPSWLKVNKIYYPD
FIVDPKKAAVWEITGAEFKSEAHTADGISIRFPRCTRIRDDKDWKSATNLPQLKELYQLSKEKADFTV
VAGDEGSSTTGGSEENKGPSGSAVSRKAPSKPSASTKKAEGKLSNSNSKDGNMQTAKPSAMKVGEKLAT
KSSPVKVGKRAAETLCQTKVLLDIFTGVRLYLPSTPDFSRLRRYFVAFDGDVLVQEFDMTSATHVLG
SRDKNPAAQQVSPewiwacIRKRLVAPC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	108.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_039269
Locus ID:	3980
UniProt ID:	P49916
RefSeq Size:	3722
Cytogenetics:	17q12
RefSeq ORF:	3027
Synonyms:	LIG2; LIG3alpha
Summary:	This gene is a member of the DNA ligase family. Each member of this family encodes a protein that catalyzes the joining of DNA ends but they each have a distinct role in DNA metabolism. The protein encoded by this gene is involved in excision repair and is located in both the mitochondria and nucleus, with translation initiation from the upstream start codon allowing for transport to the mitochondria and translation initiation from a downstream start codon allowing for transport to the nucleus. Additionally, alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome
Protein Pathways:	Base excision repair

Product images:



Coomassie blue staining of purified LIG3 protein (Cat# [TP309530]). The protein was produced from HEK293T cells transfected with LIG3 cDNA clone (Cat# [RC209530]) using MegaTran 2.0 (Cat# [TT210002]).