

Product datasheet for **TP309534L**

HSD17B7 (NM_016371) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human hydroxysteroid (17-beta) dehydrogenase 7 (HSD17B7), 1 mg
Species: Human
Expression Host: HEK293T
Expression cDNA Clone or AA Sequence: >RC209534 protein sequence
Red=Cloning site **Green**=Tags(s)

MRKVVLITGASSGIGLALCKRLLAEDDELHLCLACRNMSKAEAVCAALLASHPTAEVTIVQVDVSNLQSV
FRASKELKQRFQRLDCIYLNAGIMPNPQLNIKALFFGLFSRKVIHMFSTAEGLLTQGDKITADGLQEVFE
TNVFGHFILIRELEPLLCHSDNPSQLIWTSSRSARKSNFSLEDFQHSKGKEPYSSSKYATDLLSVALNRN
FNQQGLYSNVACPGTALTNLTYGILPPFIWTLMLPAILLRFFANAFTLTPYNGTEALVWLFHQKPESLN
PLIKYLSATTGFRNYIMTQKMDLDEDTAEKIFYQKLELEKHIRVTIQKTDNQARLSGSL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 38 kDa
Concentration: >0.05 µg/µL as determined by microplate BCA method
Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage: Store at -80°C.
Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq: [NP_057455](#)
Locus ID: 51478



[View online »](#)

UniProt ID: [P56937](#), [A0A024R913](#)

RefSeq Size: 1533

Cytogenetics: 1q23.3

RefSeq ORF: 1023


Synonyms: PRAP; SDR37C1

Summary: HSD17B7 encodes an enzyme that functions both as a 17-beta-hydroxysteroid dehydrogenase (EC 1.1.1.62) in the biosynthesis of sex steroids and as a 3-ketosteroid reductase (EC 1.1.1.270) in the biosynthesis of cholesterol (Marijanovic et al., 2003 [PubMed 12829805]).[supplied by OMIM, May 2010]

Protein Families: Druggable Genome

Protein Pathways: Androgen and estrogen metabolism, Metabolic pathways, Steroid biosynthesis

Product images:



Coomassie blue staining of purified HSD17B7 protein (Cat# [TP309534]). The protein was produced from HEK293T cells transfected with HSD17B7 cDNA clone (Cat# [RC209534]) using MegaTran 2.0 (Cat# [TT210002]).