

Product datasheet for **TP309574M**

CPNE3 (NM_003909) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human copine III (CPNE3), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209574 protein sequence Red =Cloning site Green =Tags(s)

MAAQCVTKVALNVSCANLLDKDIGSKSDPLCVLFLNTSGQQWYEVERTERIKNCLNPQFSKTFIIDYYFE
VVQKLKFGVYDIDNKTIELSDDDFLGECECTLGQIVSSKCLRPLVMKTGRPAGKGSITISAEIKNRV
VLFEMEARKLDNKDLFGKSDPYLEFHKQTSNGWLMVHRTEVVKNNLNPVWRPFKISLNSLCYGDMDKTI
KVECYDYDNDGSHDLIGTFQTTMTKLKEASRSPVEFECINEKKRQKKKSYKNSGVISVKQCEITVECTF
LDYIMGGCQLNFTVGVDFGSDPRSPDSLHYISPNGVNEYLTALWSVGLVIQDYDADKMFPAFGFGAQ
IPPQWQVSHEFPMNPNPNPCNGIQGIVEAYRSLPQIKLYGPTNFSPINHVAFAAAATQQQTASQY
FVLLIITDGVITDLDETRQAIVNASRLPMSIIIVGVGGADFSAMEFLDGDGGSLRSPLGEVAIRDIVQFV
PFRQFQNAPEALAQCVLAEIPQQVGVYFNTYKLLPPKNPATKQQKQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	60 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online >](#)

RefSeq: [NP_003900](#)

Locus ID: 8895

UniProt ID: [Q75131](#), [A0A024R994](#)

RefSeq Size: 4883


Cytogenetics: 8q21.3

RefSeq ORF: 1611

Synonyms: CPN3; PRO1071

Summary: Calcium-dependent membrane-binding proteins may regulate molecular events at the interface of the cell membrane and cytoplasm. This gene encodes a protein which contains two type II C2 domains in the amino-terminus and an A domain-like sequence in the carboxy-terminus. The A domain mediates interactions between integrins and extracellular ligands. [provided by RefSeq, Aug 2008]

Product images:



Coomassie blue staining of purified CPNE3 protein (Cat# [TP309574]). The protein was produced from HEK293T cells transfected with CPNE3 cDNA clone (Cat# [RC209574]) using MegaTran 2.0 (Cat# [TT210002]).