

## Product datasheet for TP309613

### C1orf94 (NM\_032884) Human Recombinant Protein

#### Product data:

**Product Type:** Recombinant Proteins  
**Description:** Recombinant protein of human chromosome 1 open reading frame 94 (C1orf94), transcript variant 2, 20 µg  
**Species:** Human  
**Expression Host:** HEK293T  
**Expression cDNA Clone or AA Sequence:** >RC209613 protein sequence  
**Red**=Cloning site **Green**=Tags(s)

MPVISSRQDCDSATSTVTDILCAAQVSKSGTDRGRILGDSNLQVSKLLSQFPLKSTETSKVPDNKNVL  
DKTRVTKDFLQDNLFSGPGPKPTGLSPFLLPPRPPPARPDKLPELPAQKRQLPVFAKICSKPKADPAV  
ERHHLMEWSPGKPKKGGQSLFSLQWPQSQKDACGEEGCCDAVGTASLTLPPKKPTCPAEKNLLYEFLG  
ATKNPSGQPRLRNKVEVDGPELKFNPVTVADKNNPKYTGNVFTPHFPTAMTSATLNQPLWLNLYPPPP  
VFTNHSTFLQYQGLYPQQAARMPYQQALHPQLGCYSQQVMPYNPQQMGQQIFRSSYTPLLSYIPFVQPNY  
PYPQRTPPKMSANPRDPPLMAGDGPQYLFQYGGFGSTSGGPLMHSPYFSSSGNGINF


**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

**Tag:** C-Myc/DDK  
**Predicted MW:** 44.7 kDa  
**Concentration:** >0.05 µg/µL as determined by microplate BCA method  
**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining  
**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol  
**Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.  
**Note:** For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.  
**Storage:** Store at -80°C.  
**Stability:** Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.  
**RefSeq:** [NP\\_116273](#)



[View online »](#)

Locus ID: 84970  
UniProt ID: [Q6P1W5](#)  
RefSeq Size: 2136  
Cytogenetics: 1p35.1  
RefSeq ORF: 1224

**Product images:**

Coomassie blue staining of purified C1orf94 protein (Cat# TP309613). The protein was produced from HEK293T cells transfected with C1orf94 cDNA clone (Cat# [RC209613]) using MegaTran 2.0 (Cat# [TT210002]).