

Product datasheet for **TP309616L**

ADH1C (NM_000669) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human alcohol dehydrogenase 1C (class I), gamma polypeptide (ADH1C), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209616 protein sequence Red =Cloning site Green =Tags(s)
	MSTAGKVIKCKAAVLWELKKPFSIEEVEVAPPKAHEVRIKMVAAGICRSDEHVVSGNLVTPLPVILGHEA AGIVESVGEVTTVKPGDKVIPLFTPQCGKCRICKNPESNYCLKNDLGNPRGTLQDGTRRFTCSGKPIHH FVGVSTFSQYTVVDENAVAKIDAASPLEKVCLIGCGFSTGYGSAVKVAKVTPGSTCAVFLGGVGLSVVM GCKAAGAARIIVDINKDKFAKAKELGATECINPQDYKKPIQEVLKEMTDGGVDFSEVIGQLDTMMASL LCCHEACGTSVIVGVPPDSQNLNINPMLLLTGRTWKGAIFGGFKSKESVPKLVADFMAKKFSLDALITNV LPFEKINEGFDLLRSGKSIRTVLTF
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	39.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_000660</u>



[View online »](#)


Locus ID:	126
UniProt ID:	P00326
RefSeq Size:	1769
Cytogenetics:	4q23
RefSeq ORF:	1125
Synonyms:	ADH3

Summary: This gene encodes class I alcohol dehydrogenase, gamma subunit, which is a member of the alcohol dehydrogenase family. Members of this enzyme family metabolize a wide variety of substrates, including ethanol, retinol, other aliphatic alcohols, hydroxysteroids, and lipid peroxidation products. Class I alcohol dehydrogenase, consisting of several homo- and heterodimers of alpha, beta, and gamma subunits, exhibits high activity for ethanol oxidation to acetaldehyde, thus playing a major role in ethanol catabolism. Three genes encoding alpha, beta and gamma subunits are tandemly organized in a genomic segment as a gene cluster. An association between ADH1C polymorphism and alcohol dependence has not been established. [provided by RefSeq, Sep 2019]

Protein Families: Druggable Genome

Protein Pathways: Drug metabolism - cytochrome P450, Fatty acid metabolism, Glycolysis / Gluconeogenesis, Metabolic pathways, Metabolism of xenobiotics by cytochrome P450, Retinol metabolism, Tyrosine metabolism

Product images:



Coomassie blue staining of purified ADH1C protein (Cat# [TP309616]). The protein was produced from HEK293T cells transfected with ADH1C cDNA clone (Cat# [RC209616]) using MegaTran 2.0 (Cat# [TT210002]).