

Product datasheet for **TP309650M**

OTUB2 (NM_023112) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human OTU domain, ubiquitin aldehyde binding 2 (OTUB2), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209650 protein sequence Red =Cloning site Green =Tags(s)

MSETSFNLISEKCDILSILRDHPENRIYRRKIEELSKRFTAIRKTKGDGNCFYRALGYSYLESLLGKSRE
IFKFKERVLQTPNDLLAAGFEEHKFRNFFNAFYVWELVEKDGVSLLKVFNDQSASDHIVQFLRLLTS
AFIRNRADFFRHFIDEEMDIKDFCTHEVEPMATECDHIQITALSQALSIALQVEYVDEMDTALNHHVFPE
AATPSVYLLYKTSHYNILYAADKH

SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	27 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_075601</u>
Locus ID:	78990
UniProt ID:	<u>Q96DC9</u>



[View online »](#)

RefSeq Size: 3873

Cytogenetics: 14q32.12


RefSeq ORF: 702

Synonyms: C14orf137; OTB2; OTU2

Summary: This gene encodes one of several deubiquitylating enzymes. Ubiquitin modification of proteins is needed for their stability and function; to reverse the process, deubiquitylating enzymes remove ubiquitin. This protein contains an OTU domain and binds Ubal (ubiquitin aldehyde); an active cysteine protease site is present in the OTU domain. [provided by RefSeq, Aug 2011]

Protein Families: Protease

Product images:



Coomassie blue staining of purified OTUB2 protein (Cat# [TP309650]). The protein was produced from HEK293T cells transfected with OTUB2 cDNA clone (Cat# [RC209650]) using MegaTran 2.0 (Cat# [TT210002]).