

Product datasheet for TP309712

ZNF524 (NM_153219) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human zinc finger protein 524 (ZNF524), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209712 protein sequence Red =Cloning site Green =Tags(s) MDTPSPDPLPSPLPGEEEEKPLALSPPVPRGRRRPGGATSSNRTLKASLHRKRGRPPKSGQEPPLVQVQ GVTAPVGSSGGSDLLIDQGVPTVSEGSAAAGPEGSGPRKAPHFCPVCLRAFPYLSDLERHSISHSELK PHQCKVCGKTFKRSSHLLRRHCNIHAGLRPFRCPLCTRRFREAGELAHHRVHSGERPYQCPICRLRFTEA NTLRRHAKRKHPEAMGVPLCAPDPGSEPPWDEEGIPATAGAEETEGKGPEA TR TRPLEQKLISEEDLA NDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	28.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_694951</u>
Locus ID:	147807
UniProt ID:	<u>Q96C55</u>


[View online »](#)


RefSeq Size: 1135

Cytogenetics: 19q13.42

RefSeq ORF: 792

Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified ZNF524 protein (Cat# TP309712). The protein was produced from HEK293T cells transfected with ZNF524 cDNA clone (Cat# [RC209712]) using MegaTran 2.0 (Cat# [TT210002]).