

Product datasheet for TP309780L

PUS7L (NM_031292) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human pseudouridylate synthase 7 homolog (S. cerevisiae)-like (PUS7L), transcript variant 3, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209780 protein sequence Red =Cloning site Green =Tags(s)
	MEEDTDYRIRFSSLCFFNDHVGFGHTIKSSPSDFVIEIDEQGQLVNKTIDEPIFKISEIQLEPNNFPKK PKLDLQNLSELDGRNQEVHTLIKYTGDQNHQSGSEKEDTIVDGTSCCEKADVLSSFLDEKTHELLNNF ACDVREKWLSTELIGLPPEFSIGRIVDKNQASLHSAIRQKFPFLVTVGKNSEIVVKPNLEYKELCHLV SEEEAFDFFKYLDAKKENSKFTFKPDTNKHAKAVHHFVNKKFGNLVETKSFSKMNCSSAGNPVNVTVRF REKAHKRGKRPLSECQEGKVIYTAFTLRKENLEMFEAIGFLAIKLGVIPSDFSYAGLKDKKAITYQAMWV RKVTPERLKNIEKEIEKKRMNVFNIRSVDDSLRLGQLKGNHFDIVIRNLKKQINDSANLRERIMEAIENV KKKGFVNYYGPQRFQGRKVHTDQIGLALLKNEMMKAIKFLTPEDLDDPVNRKKYFLQTEDAKGTLVL MPEFKVRERALLEALHRFGMTEEGCIQAWFSLPHSMRIFYVHAYTSKIWNEAVSYRLETYGARVVQGDV CLDEDIDDENFQNSKIHLVTEEEGSANMYAIHQVLPVLGYNIQYPKNKVGQWYHDILSRDGLQTCRFKV PTLKLNIPGCRYQLKHPCNLSYQLMEDHDIDVKTGSHIDETALLISFDLDASCYATVCLKEIMKHD V TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	80.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



View online »

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_112582
Locus ID:	83448
UniProt ID:	Q9H0K6
RefSeq Size:	3985
Cytogenetics:	12q12
RefSeq ORF:	2103
Summary:	Pseudouridylate synthase that catalyzes pseudouridylation of RNAs.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified PUS7L protein (Cat# [TP309780]). The protein was produced from HEK293T cells transfected with PUS7L cDNA clone (Cat# [RC209780]) using MegaTran 2.0 (Cat# [TT210002]).