

Product datasheet for **TP309811M**

FAM156A (NM_014138) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human family with sequence similarity 156, member A (FAM156A), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209811 protein sequence Red =Cloning site Green =Tags(s)
	 MDPLQKRNPASPSKSSPMTAAETSQEGPAPSQPSYSEQPMMGLSNLSPGPGHSQAVPLPEGLLHQRYREE KTLEERRWERLEFLQRKKAFLRHVRRRHRDHMAPYAVGREARISPLGDRSQNRFQCECRYCQSHRPNLSG IPGESNRAPHPSSWETLVQGLSGLTSLGNTNQPGLPEAALQPQETEEKRQRERQQESKIMFQRLKQWL EEN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	24.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_054857</u>
Locus ID:	29057



[View online »](#)

UniProt ID: [Q8NDB6](#)
RefSeq Size: 2248
Cytogenetics: Xp11.22
RefSeq ORF: 639
Synonyms: PRO0659; TMEM29

Product images:

Coomassie blue staining of purified FAM156A protein (Cat# [TP309811]). The protein was produced from HEK293T cells transfected with FAM156A cDNA clone (Cat# [RC209811]) using MegaTran 2.0 (Cat# [TT210002]).