

Product datasheet for TP309816

DNA Polymerase lambda (POLL) (NM_013274) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human polymerase (DNA directed), lambda (POLL), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209816 protein sequence Red =Cloning site Green =Tags(s)

MDPRGILKAFPKRQKIHADASSKVLAKIPRREEGEEAEWLSSLRAHVVRTGIGRARAELFEKQIVQHGG
QLCPAQGPVTHIVDEGM DYERALLRLPQLPPGAQLVKS AWLSLCLQERRLVDVAGFSIFIPSRYLD
HPQPSKAEQDASIPPGTHEALLQTALSPPPPTRPVSPQKAKEAPNTQAQPISDDEASDGEETQVSAAD
LEALISGHYPTSLEGDCEPSPAPAVLDKWVCAQPSSQKATNHNHLHITEKLEVLAKAYSVQGDKWRALGYA
KAINALKSFHKPVTSYQGAC SIPGIGKRMAEKIIELESGLRKL DHISESVPVLELFSNIW GAGTKTAQ
MWYQQGFRSLEDIRSQASLTTQQAIGLKHYSDFLERMPREEATEIEQTVQKAAQAFNSGLLCVACGSYRR
GKATCGDVDVLITHPDGRSHRGIFSRLLDSL RQEGFLTDDLVSQEENGQQQKYLGVCR LPPGPRRHRRLD
IIVPYSEFACALLYFTGSAHFNRSMRALAKTKGMSLSEHALSTAVVRNTHGCKVGPGRVLP TPEKDFV
RLLGLPYREPAERDW

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	63.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_037406
Locus ID:	27343
UniProt ID:	Q9UGP5
RefSeq Size:	2741
Cytogenetics:	10q24.32
RefSeq ORF:	1725
Synonyms:	BETAN; POLKAPPA
Summary:	This gene encodes a DNA polymerase. DNA polymerases catalyze DNA-template-directed extension of the 3'-end of a DNA strand. This particular polymerase, which is a member of the X family of DNA polymerases, likely plays a role in non-homologous end joining and other DNA repair processes. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2010]
Protein Families:	Druggable Genome
Protein Pathways:	Base excision repair, Non-homologous end-joining

Product images:



Coomassie blue staining of purified POLL protein (Cat# TP309816). The protein was produced from HEK293T cells transfected with POLL cDNA clone (Cat# [RC209816]) using MegaTran 2.0 (Cat# [TT210002]).