

Product datasheet for **TP309817M**

MKLN1 (NM_013255) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human muskelin 1, intracellular mediator containing kelch motifs (MKLN1), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209817 protein sequence Red =Cloning site Green =Tags(s) MAAGGAVAAAP ECRLLPYALHKWSSFSSTYLPENILVDKPNDQSSRWSSSESNYPQYLILKLERPAIVQN ITFGKYEKTHVCNLKKFKVFGGMNEENMTELLSSGLKNDYNKETFTLKHKIDEQMFPCRFIKIVPLLSWG PSFNFSIWYVELSGIDDPDIVQPCLNWYSKYREQEAIRLCLKHFRQHNYTEAFESLQKKTIALEHPMLT DIHDKLV LKGFDFACEELIEKAVNDGLFNQYISQQEYKPRWSQIIPKSTKGDGEDNRP GMRGGHQMVIDV QTETVYLFGGWDGTQDLADFWAYSVKENQWTCISRDEKENGPSARSCHKMCIDIQRRQIYTLGRYLDSS VRNSKSLKSDFYRYDIDTNTWMLLSEDTAADGGPKLVFDHQMCM DSEKHMIYTFGGRI LTCNGSVDDSR A SEPQFSGLF AFNCQCQ TWKLLREDSCNAGPEDIQSRIGHCMLFHSKNRCLYVFGGQRSKTYLNDFFSYDV DSDHVDIISDGTKKDSGMVPMTGFTQRTIDPELNEIHVLSGLSKDKEKREENVRNSFWIYDIVRNSWSC VYKNDQAAKDNPTKSLQEEPCPRFAHQLVYDELHKVHYLFGGNPGKSCSPKMRLDDFWSLKLCRPSKD Y LLRHCKYLIRKHRFEEKAQVDPLSALKYLQNDLYITVDHSDPEETKEFQLLASALFKSGSDF TALGFSDV DHTYAQRTQLFDTLVNFFPD SMTPPKGNLVDLITL TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	84.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.


[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_037387
Locus ID:	4289
UniProt ID:	Q9UL63
RefSeq Size:	11194
Cytogenetics:	7q32.3
RefSeq ORF:	2205
Synonyms:	TWA2
Summary:	Muskelin is an intracellular protein that acts as a mediator of cell spreading and cytoskeletal responses to the extracellular matrix component thrombospondin I (MIM 188060) (Adams et al., 1998 [PubMed 9724633]).[supplied by OMIM, Mar 2008]

Product images:



Coomassie blue staining of purified MKLN1 protein (Cat# [TP309817]). The protein was produced from HEK293T cells transfected with MKLN1 cDNA clone (Cat# [RC209817]) using MegaTran 2.0 (Cat# [TT210002]).