

Product datasheet for **TP309854L**

VPS39 (NM_015289) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human vacuolar protein sorting 39 homolog (S. cerevisiae) (VPS39), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209854 representing NM_015289 Red =Cloning site Green =Tags(s)

MHDAFEPVPILEKLPLQIDCLAAWEEWLLVGTKQGHLLYRIRKDVGCNRFEVTLEKSNKNFSKKIQQIH
 VVSQFKILVSLLENNIYVHDLTFQQITTVSKAKGASLFTCDLQHTETGEEVLRMCVAVKKKLQLYFWKD
 REFHELQGDVSPDVPKSMAWCENSICVGFKRDYYLIRVDGKGSIKELFPTGKQLEPLVAPLADGKVAVG
 QDDLTVVLNEEGICTQKCALNWDIPVAMEHQPPYIIAVLPRYVEIRTFEPRLLVQSIELQRPRFITSGG
 SNIIYVASNHFWRLIPVPMATQIQQLQDKQFELALQLAEMKDDSDSEKQQQIHHIKNLYAFNLFQCKR
 FDESMQVFAKLGTDPTHVMGLYPDLLPTDYRKQLQYPNPLPVLSGAELEKAHLALIDYLTQKRSQVLVKKL
 NDS DHQSSTPLMEGTPTIKSKKKLLQIIDTLLKCYLHTNVALVAPLLRLENNHCHIEESEHVLKKAHK
 YSELIILEYKKGLEKALQVLVDQSKKANSPLKGHERTVQYLQHLGTENLHLIFSYSVWVLRDFPEDGLK
 IFTEDLPEVESLPRDRVLFGLIENFKGLAIPYLEHIIHVWEETGSRFHNCLIQLYCEKVQGLMKEYLLSF
 PAGKTPVPAGEEEGELGEYRQKLLMFLEISSYDPRGLICDFPFDGLLEERALLGRMGKHEQALFIYVH
 ILKDTRMAEEYCHKHYDRNKDGNKDVYLSLLRMYLSPPSIHCLGPIKLELLEPKANLQAALQVLELHHSK
 LDTTKALNLLPANTQINDIRIFLEKVLEENAKKRFNQVLKNLLHAEFLRVQEERILHQVVKCIITEEKV
 CMVCKKKIGNSAFARYPNGVWVHYFCSKEVNPADT


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	100.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.


[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_056104
Locus ID:	23339
UniProt ID:	Q96JC1
RefSeq Size:	4851
Cytogenetics:	15q15.1
RefSeq ORF:	2625
Synonyms:	hVam6p; TLP; VAM6
Summary:	This gene encodes a protein that may promote clustering and fusion of late endosomes and lysosomes. The protein may also act as an adaptor protein that modulates the transforming growth factor-beta response by coupling the transforming growth factor-beta receptor complex to the Smad pathway. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]
Protein Families:	Druggable Genome

Product images:



Coomassie blue staining of purified VPS39 protein (Cat# [TP309854]). The protein was produced from HEK293T cells transfected with VPS39 cDNA clone (Cat# [RC209854]) using MegaTran 2.0 (Cat# [TT210002]).