

Product datasheet for **TP309858M**

C22orf40 (CDPF1) (NM_207327) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 22 open reading frame 40 (C22orf40), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209858 protein sequence Red =Cloning site Green =Tags(s)
	MASHVECRPLGVFECELCTLTAPYSYVGQKPPNTQSMVLLEESYVMKDPFTSDKDRFLVLGSCCSLCSRL VCVGPECSLFYSKRFLPCVRENINAFQEIQRQDLEKRRKAPSKRTPSQPGSRT
	TR TRPLE QKLISEEDLA ANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	13.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_997210
Locus ID:	150383
UniProt ID:	Q6NVV7
RefSeq Size:	1502
Cytogenetics:	22q13.31




[View online »](#)

RefSeq ORF: 369

Synonyms: C22orf40

Product images:



Coomassie blue staining of purified CDPF1 protein (Cat# [TP309858]). The protein was produced from HEK293T cells transfected with CDPF1 cDNA clone (Cat# [RC209858]) using MegaTran 2.0 (Cat# [TT210002]).