

Product datasheet for TP309863

PWWP domain containing protein 2B (PWWP2B) (NM_138499) Human Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant protein of human PWWP domain containing 2B (PWWP2B), transcript variant 1, 20 µg |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | >RC209863 protein sequence Red=Cloning site Green=Tags(s) |

MQLGSSSPPPACGVQPPETTGPPEPPPLVPPLPAGSLPPYPPYFEGAPFPHPLWLRDITYKLWVPQPPRT
 IKRTRRRLSRNRDPGRLILSTIRLRPRQVLCEKCKSTLSPPEASPGPPAAPRARRRLGSGPDRELKPEE
 PENGEPTAAATARRSKRERREEDRAPAEQVPRSPVIKISYSTPQGKGEVVKIPSRVHGSLEPFRPQQAPQ
 DDGSQDPEVLDRESRDRPSCAPSASIPKLLTRPVPAGADLPPPKIRLKPRLGDSEHEPVYRAELVGEL
 NGYLRDSSPAPCADGPAGGLADLSSGSSGEDDDFKSCPQGPQGREGLAFLVSCPEGRADCASESACSSDS
 LDEARSSGSEGTADTGDLSPGHGASAPSVSREARQTVPPLTVRLHTQVSSECITEDGRTVAVGDIWVGK
 IHGFPWWPARVLDISLGQKEDGEPWSWREAKVSWFGSPPTTSFLSISKLSPFSEFFKLRFNRRKKKGMRYKAI
 TEAANAARHVAPEIRELLTQFET


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

| | |
|----------------|--|
| Tag: | C-Myc/DDK |
| Predicted MW: | 63.8 kDa |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol |
| Preparation: | Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps. |
| Note: | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage: | Store at -80°C. |



[View online »](#)

| | |
|---------------|---|
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| RefSeq: | NP_612508 |
| Locus ID: | 170394 |
| UniProt ID: | Q6NUJ5 |
| RefSeq Size: | 2626 |
| Cytogenetics: | 10q26.3 |
| RefSeq ORF: | 1539 |
| Synonyms: | bA432J24.1; pp8607; PWWP2 |

Product images:

Coomassie blue staining of purified PWWP2B protein (Cat# TP309863). The protein was produced from HEK293T cells transfected with PWWP2B cDNA clone (Cat# [RC209863]) using MegaTran 2.0 (Cat# [TT210002]).