

Product datasheet for TP309879L


SPLUNC3 (BPIFA3) (NM_001042439) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 20 open reading frame 71 (C20orf71), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209879 protein sequence Red=Cloning site Green=Tags(s)
	<p>MMCPLWRLIFLGLLALPLAPHKQPWPGLAQAHARDNKSTLARIIAQGLIKHNAESRIQNIHFGDRLNASA QVAPGLVGWLISGRKHQQQESRSFDNNIVKMCAMHSIVVEFWLEKDEFGRDLVIGKCDAPSSVHVAI LTEAIPPKMNQFLYNLKENLQKVLPHMVESQVCPLIGEILGQLDVKLLKSLIEQAAHEPTHHETSQPSA CQAGESPS</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	24.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001035904</u>
Locus ID:	128861


[View online »](#)

UniProt ID: Q9BQP9
 RefSeq Size: 1087
 Cytogenetics: 20q11.21
 RefSeq ORF: 654
 Synonyms: C20orf71; SPLUNC3

Product images:


Coomassie blue staining of purified BPIFA3 protein (Cat# [TP309879]). The protein was produced from HEK293T cells transfected with BPIFA3 cDNA clone (Cat# [RC209879]) using MegaTran 2.0 (Cat# [TT210002]).