

Product datasheet for TP309891L

MARCHF8 (NM_145021) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human membrane-associated ring finger (C3HC4) 8 (MARCH8), transcript variant 8, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209891 protein sequence Red =Cloning site Green =Tags(s) MSMP LHQISAIP SQDAISARVYRSKTK EKEREEQNEKTLGHFMSHSSNISKAGSPPSASAPAVSSFSRT SITPSSQDICRICHCEGDDESPLITPCHCTGSLHFVHQACLQQWIKSSDTRCCELCKYEFIMETKLKPLR KWEKLQMTSSERRKIMCSVTFHVIAITCVVWSLYVLIDRTAEIKQGQATGILEWPFWTKLVVVAIGFTG GLLFMVYQCKVYVQLWKRLKAYNRVIYVQNC PETS KKNIFEKSPLTEPNFENKHGHGICHSDTNSSCCTE PEDTGAEIIHV TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	32.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_659458</u>
Locus ID:	220972


[View online »](#)

UniProt ID:	<u>Q5T0T0</u>
RefSeq Size:	5300
Cytogenetics:	10q11.21-q11.22
RefSeq ORF:	873
Synonyms:	c-MIR; CMIR; MARCH-VIII; MARCH8; MIR; RNF178
Summary:	MARCH8 is a member of the MARCH family of membrane-bound E3 ubiquitin ligases (EC 6.3.2.19). MARCH enzymes add ubiquitin (see MIM 191339) to target lysines in substrate proteins, thereby signaling their vesicular transport between membrane compartments. MARCH8 induces the internalization of several membrane glycoproteins (Goto et al., 2003 [PubMed 12582153]; Bartee et al., 2004 [PubMed 14722266]).[supplied by OMIM, Apr 2010]
Protein Families:	Druggable Genome, Transmembrane

Product images:



Coomassie blue staining of purified MARCH8 protein (Cat# [TP309891]). The protein was produced from HEK293T cells transfected with MARCH8 cDNA clone (Cat# [RC209891]) using MegaTran 2.0 (Cat# [TT210002]).