

Product datasheet for TP309893

FLJ43980 (NM_001004299) Human Recombinant Protein

Product data:


| | |
|---------------------------------------|---|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant protein of human FLJ43980 protein (FLJ43980), 20 µg |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | >RC209893 representing NM_001004299 Red =Cloning site Green =Tags(s) |
| | MRRSSGFGGQKKGQGPSCSFTGCWCCRGDDVAESDDSPFAQCGYNIQEKHLGKLHRAASRGEVSKVECILS SGSADLDERDKKKRTALHLACANGHPEVALLVDRGCQLDVFNDNKNRTALLKAVQCQEEECATILLEHGA DPDLPDVYGNLTLHYAIYNEDIPMTKLLHHANIESANKDELTPFLLAVHEQKQMEDFLRKQKENLTA VKLESIHQVMSEYKENETPRNPQNSNPEGTSNKMACLGEGAAGAKVDEIPGNPVTRLFNKPSIDDSRPMS ANEDFDFDTEEKATEPANGKRQNGMGIIESAPQEHTNNENI |
| | TRTRPLEQKLISEEDLAANDILDYKDDDDKV |
| Tag: | C-Myc/DDK |
| Predicted MW: | 35.9 |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol |
| Preparation: | Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps. |
| Note: | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage: | Store at -80°C. |
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| RefSeq: | <u>NP_001004299</u> |
| Locus ID: | 124149 |



[View online »](#)

RefSeq Size: 3722
Cytogenetics: 16q11.2
RefSeq ORF: 963
Synonyms: MGC87661, DKFZp781D1722

Product images:



Coomassie blue staining of purified FLJ43980 protein (Cat# TP309893). The protein was produced from HEK293T cells transfected with FLJ43980 cDNA clone (Cat# [RC209893]) using MegaTran 2.0 (Cat# [TT210002]).