

Product datasheet for TP309909L

14-3-3 zeta (YWHAZ) (NM_003406) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, zeta polypeptide (YWHAZ), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209909 representing NM_003406 Red=Cloning site Green=Tags(s)

MDKNELVQKAKLAEQAERYDDMAACMKSVTEQGAELSNEERNLLSVAYKNVVGARRSSWRVSSIEQKTE
GAEKKQQMAREYREKIETELRDICNDVLSLLEKFLIPNASQAESKVFLKMKGDYYRYLAEVAAGDDKKG
IVDQSQQAYQEAFEISKEMQPHTPIRLGLALNFSVFYIEILNSPEKACSLAKTAFDEAIAELDTLSEES
YKDSTLIMQLLRDNLTLWTSDTQGDEAEAGEGGEN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	27.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_003397
Locus ID:	7534



[View online »](#)

UniProt ID: [P63104](#), [D0PNI1](#)

RefSeq Size: 2834

Cytogenetics: 8q22.3


RefSeq ORF: 735

Synonyms: 14-3-3-zeta; HEL-S-3; HEL-S-93; HEL4; KCIP-1; POPCHAS; YWHAD

Summary: This gene product belongs to the 14-3-3 family of proteins which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 99% identical to the mouse, rat and sheep orthologs. The encoded protein interacts with IRS1 protein, suggesting a role in regulating insulin sensitivity. Several transcript variants that differ in the 5' UTR but that encode the same protein have been identified for this gene. [provided by RefSeq, Oct 2008]

Protein Pathways: Cell cycle, Neurotrophin signaling pathway, Oocyte meiosis, Pathogenic Escherichia coli infection

Product images:



Coomassie blue staining of purified YWHAZ protein (Cat# [TP309909]). The protein was produced from HEK293T cells transfected with YWHAZ cDNA clone (Cat# [RC209909]) using MegaTran 2.0 (Cat# [TT210002]).