

Product datasheet for **TP309916M**

RETREG1 (NM_019000) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human family with sequence similarity 134, member B (FAM134B), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209916 protein sequence Red =Cloning site Green =Tags(s)
	<p>MPEGEDFGPGKSWEVINSKPDERPRLSHCIAESWMNFSIFLQEMSLFKQQSPGKFLLVCSVCTFFFTILG SYIPGVILSYLLLLCAFLCPLFKCNDIGQKIYSKIKSVLLKDFGIGEYINQKKRERSEADKEKSHKDDSD ELDFSALCPKISLTVAAKELSVSDTDVSEVSWTDNGTFNLSEGYTPQTDTSDDLDRPSEEVFSRDLSDFP SLENGMGTNDEDELSLGLPTELKRKKEQLDSGHRPSKETQSAAGLTLPLNSDQTFHLMNSNLAGDVITAAV TAAIKDQLEGVQQALSQAAPIPEEDTDTEEGDDFELLDQSELDQIESELGLTQDQEAQAQNKSSGFLS NLLGGH</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	39.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_061873</u>




[View online »](#)

Locus ID: 54463
UniProt ID: [Q9H6L5](#)
RefSeq Size: 3184
Cytogenetics: 5p15.1
RefSeq ORF: 1068
Synonyms: FAM134B; JK-1; JK1

Summary: The protein encoded by this gene is a cis-Golgi transmembrane protein that may be necessary for the long-term survival of nociceptive and autonomic ganglion neurons. Mutations in this gene are a cause of hereditary sensory and autonomic neuropathy type IIB (HSAN IIB), and this gene may also play a role in susceptibility to vascular dementia. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Nov 2011]

Protein Families: Transmembrane

Product images:



Coomassie blue staining of purified RETREG1 protein (Cat# [TP309916]). The protein was produced from HEK293T cells transfected with RETREG1 cDNA clone (Cat# [RC209916]) using MegaTran 2.0 (Cat# [TT210002]).