

Product datasheet for **TP309959M**

CENPP (NM_001012267) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human centromere protein P (CENPP), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209959 protein sequence Red =Cloning site Green =Tags(s)

MDAELAEVRALQAEIAALRRACEDPPAPWEEKSRVQKSFQAIHQFNLEGWKSSKDLKNQLGHLESELSFL
STLTGINIRNHSKQTEDLTSTEMTEKSIRKVLQRHRLSGNCHMVTFQLEFQILEIQNKERLSSAVTDLNI
IMEPTCESELSEFVSRAEERKDLFMFFRSLHFFVEWFEYRKRTEFKHLKEKYPDAVYLSEGPSSCSMGIRS
ASRPGFELVIVWRIQIDEDGKVFPKLDLLTKVPQRALELDKNRAIETAPLSFRTLVLGILLGIEAALESLIK
SLCAEENN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	33 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001012267
Locus ID:	401541



[View online »](#)

UniProt ID: [Q6IPU0](#)

RefSeq Size: 3413


Cytogenetics: 9q22.31

RefSeq ORF: 864

Synonyms: CENP-P

Summary: CENPP is a subunit of a CENPH (MIM 605607)-CENPI (MIM 300065)-associated centromeric complex that targets CENPA (MIM 117139) to centromeres and is required for proper kinetochore function and mitotic progression (Okada et al., 2006 [PubMed 16622420]). [supplied by OMIM, Mar 2008]

Product images:



Coomassie blue staining of purified CENPP protein (Cat# [TP309959]). The protein was produced from HEK293T cells transfected with CENPP cDNA clone (Cat# [RC209959]) using MegaTran 2.0 (Cat# [TT210002]).