

Product datasheet for **TP309974L**

AK2 (NM_001625) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human adenylate kinase 2 (AK2), transcript variant AK2A, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC209974 protein sequence Red =Cloning site Green =Tags(s)
	<p>MAPSVPAAEPEYPKGIRAVLLGPPGAGKGTQAPRLAENFCVCHLATGDMLRAMVASGSELGKKLKATMDA GKLVSDEMVELIEKNLETPLCCKNGFLLDGFPRTVRQAEMLDDLMEKRKEKLDVIEFSIPDSSLIRIT GRLIHPKSGRSYHEEFNPPKEPMKDDITGEPLIRRSDDNEKALKIRLQAYHTQTTPLEYYRKRGIHSAI DASQTPDVVFASILAAFSKATCKDLVMFI</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	26.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001616
Locus ID:	204
UniProt ID:	P54819 , A0A140VK93



[View online »](#)

RefSeq Size: 2759

Cytogenetics: 1p35.1

RefSeq ORF: 717


Synonyms: ADK2

Summary: Adenylate kinases are involved in regulating the adenine nucleotide composition within a cell by catalyzing the reversible transfer of phosphate groups among adenine nucleotides. Three isozymes of adenylate kinase, namely 1, 2, and 3, have been identified in vertebrates; this gene encodes isozyme 2. Expression of these isozymes is tissue-specific and developmentally regulated. Isozyme 2 is localized in the mitochondrial intermembrane space and may play a role in apoptosis. Mutations in this gene are the cause of reticular dysgenesis. Alternate splicing results in multiple transcript variants. Pseudogenes of this gene are found on chromosomes 1 and 2.[provided by RefSeq, Nov 2010]

Protein Families: Druggable Genome

Protein Pathways: Metabolic pathways, Purine metabolism

Product images:



Coomassie blue staining of purified AK2 protein (Cat# [TP309974]). The protein was produced from HEK293T cells transfected with AK2 cDNA clone (Cat# [RC209974]) using MegaTran 2.0 (Cat# [TT210002]).