

## Product datasheet for **TP31009M**

### HGD (NM\_000187) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human homogentisate 1,2-dioxygenase (homogentisate oxidase) (HGD), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210009 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MAELKYISGFGNECSSEDPKCPGSLPEGQNNPQVCPYNLYAEQLSGSAFTCPRSTNKRSLYRILPSVSH  
KPFESIDEGHVTHNWDEVDPDPNQLRWKPFEPKASQKKVDFVSLHTLCGAGDIKSNNGLAIHIFLCNT  
SMENRCFYNSDGDFLVLPQKGNLLIYTEFGKMLVQPNEICVIQRGMRFSDVFEETRGYILEVYGVHFL  
PDLGPIGANLANPRDFLIPIAWYEDRQVPGGYTVINKYQGKLFAAKQDVSPFNVAWHGNYTPYKYNLK  
NFMVINSVAFDHADPSIFTVLTAKSVRPGVAIADFVIFPPRWGVADKTFRPPYYHRNCMSEFMGLIRGHY  
EAKQGGFLPGGGSLHSTMPHGPDADCFEKASKVKLAPERIADGTMAFMFESSLSLAVTKWGLKASRCLD  
ENYHKCWEPLKSHFTPNRNPAPEN

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	49.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online >](#)

RefSeq: [NP\\_000178](#)

Locus ID: 3081

UniProt ID: [Q93099](#)

RefSeq Size: 2012

Cytogenetics: 3q13.33

RefSeq ORF: 1335


Synonyms: AKU; HGO

**Summary:** This gene encodes the enzyme homogentisate 1,2 dioxygenase. This enzyme is involved in the catabolism of the amino acids tyrosine and phenylalanine. Mutations in this gene are the cause of the autosomal recessive metabolism disorder alkaptonuria.[provided by RefSeq, May 2010]

**Protein Families:** Druggable Genome

**Protein Pathways:** Metabolic pathways, Tyrosine metabolism

### Product images:



Coomassie blue staining of purified HGD protein (Cat# [TP310009]). The protein was produced from HEK293T cells transfected with HGD cDNA clone (Cat# [RC210009]) using MegaTran 2.0 (Cat# [TT210002]).