

Product datasheet for **TP310115M**

IFNA21 (NM_002175) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human interferon, alpha 21 (IFNA21), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210115 protein sequence Red =Cloning site Green =Tags(s)
	 MALSFLLMAVLVLSYKSICSLGCDLPQTHSLGNRRALILLAQMGRISPFSCCLKDRHDFGFPQEFDGNQ FQKAQAISVLHEMIQQTENLSTKSSATWEQSLLLEKFSTELNQQLNLEACVIQEVGVEETPLMNVDSI LAVKKYFQRITLYLTEKKYSPCAWEVVRAEIMRSFSLSKIFQERLRRKE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	21.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_002166
Locus ID:	3452
UniProt ID:	P01568
RefSeq Size:	1024



[View online »](#)

Cytogenetics:	9p21.3
RefSeq ORF:	567
Synonyms:	IFN-alpha1; IeIF-F; IeIF F
Summary:	This gene is a member of the alpha interferon gene cluster on the short arm of chromosome 9. Interferons are cytokines produced in response to viral infection that mediate the immune response and interfere with viral replication. The encoded protein is a type I interferon and may play a specific role in the antiviral response to rubella virus. [provided by RefSeq, Sep 2011]
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Antigen processing and presentation, Autoimmune thyroid disease, Cytokine-cytokine receptor interaction, Cytosolic DNA-sensing pathway, Jak-STAT signaling pathway, Natural killer cell mediated cytotoxicity, Regulation of autophagy, RIG-I-like receptor signaling pathway, Toll-like receptor signaling pathway

Product images:



Coomassie blue staining of purified IFNA21 protein (Cat# [TP310115]). The protein was produced from HEK293T cells transfected with IFNA21 cDNA clone (Cat# [RC210115]) using MegaTran 2.0 (Cat# [TT210002]).