

Product datasheet for TP310308

PI 3 Kinase regulatory subunit 4 (PIK3R4) (NM_014602) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human phosphoinositide-3-kinase, regulatory subunit 4 (PIK3R4), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210308 protein sequence Red=Cloning site Green=Tags(s)

MGNQLAGIAPSQILSVESYFSDIHDFFEYDKSLGSTRFFKVARAKHREGLVVKVFAIQDPTLPLTSYKQE
LEELKIRLNSAQNCLPFQKASEKASEKAAMLFRQYVRDNLDRISTRPFLNNEIKRWIAFQILTAVDQAH
KSGVRHGDIKTENVMVTSWNWVLLTDFASFKPTYPEDNPADFNFFDTSRRRTCYIAPERFVDGGMFAT
ELEYMRDPSTPLVDLNSNQRTGELKRAMDIFSAGCVIAELFTEGVPLFDLSQLLAYRNGHFFPEQVLNK
IEDHSIRELVTQMIHREPDKRLEAEDYLKQQRGNAPPEIFYTFLQPYMAQFAKETFLSADERILVIRKDL
GNIIHNLCCGHDLPEKAEGEPKENGLVILVSVITSCLQTLKYCDSKLALELILHLAPRLSVEILLDRITP
YLLHFSNDSVPRVRAEALRTLTKVLALVKEVPRNDINIYPEYILPGIAHLAQDDATIVRLAYAENIALLA
ETALRFLELVQLKLNLMENDPNNEEIDEVTHPNGNYDELQALHEMVQKQVVTLLSDPENIVKQTLMENG
ITRLCVFFGRQKANDVLLSHMITFLNDKNDWHLRGAFFDSIVGVAAYVGVQSSSILKPLLQQGLSDAEFF
VIVKALYALTCMCQLGLLQKPHVYEFASDIAPFLCHPNLWIRYGAVGFITVVARQISTADVCKLMPYLD
PYITQPIIQUIERKLVLLSVLKEPVSRIFDYALRSKDITSLFRHLHMRQKKRNGSLPDCPPPEDPAIAQL
LKKLLSQGMTEEEEDKLLALKDFMMSKNKAKANIVDQSHLHDSSQKGVIDLALGITGRQVDLVKTKQEP
DDKRARKHVQDSNVNEEWKSMFGSLDPPNMPQALPKGSDQEVITGKPPRSESSAGICVPLSTSSQVPE
VTTVQNKKPVIPVLSSTILPSTYQIRITTCKTELQQLIQKREQCNAERIAKQMMENAEWESKPPPPGWR
PKGLLVAHLHEHKS AVNRIRVSDEHSLFATCSNDGTVKIWNWSQKMEGKTTTTRSILTYSRIGGRVKTLLF
CQGSHYLAIASDNGAVQLLGIEASKLPKSPKIHPLQSRILDQKEDGCVVDMHMFNSGAQSVLAYATVNGS
LVGWDLRSSNAWTLKHDLKSLITSAVDIHCWLCIGTSSGTMACWDMRFQLPISSHCHPSRARIRRL
SMHPLYQSWVIAAVQGNNEVSMWDMETGDRRFTLWASSAPLSELQSPHVSVHGICYSPADGNPILLTAG
SDMKIRFWDLAYPERSYVAGSTSSPSVSYRKKIIEGTEVVQEIQNKQKQVGPSDDTPRRGPESLPVGHHD
IITDVATFQTTQGFIVTASRDGIVKVVK


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	152.9 kDa



[View online »](#)

Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_055417
Locus ID:	30849
UniProt ID:	Q99570
RefSeq Size:	5078
Cytogenetics:	3q22.1
RefSeq ORF:	4074
Synonyms:	p150; VPS15
Summary:	Regulatory subunit of the PI3K complex that mediates formation of phosphatidylinositol 3-phosphate; different complex forms are believed to play a role in multiple membrane trafficking pathways: PI3KC3-C1 is involved in initiation of autophagosomes and PI3KC3-C2 in maturation of autophagosomes and endocytosis. Involved in regulation of degradative endocytic trafficking and cytokinesis, probably in the context of PI3KC3-C2 (PubMed:20643123).[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Regulation of autophagy

Product images:

Coomassie blue staining of purified PIK3R4 protein (Cat# TP310308). The protein was produced from HEK293T cells transfected with PIK3R4 cDNA clone (Cat# [RC210308]) using MegaTran 2.0 (Cat# [TT210002]).