

Product datasheet for **TP310351L**

CCDC157 (NM_001017437) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coiled-coil domain containing 157 (CCDC157), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210351 protein sequence Red =Cloning site Green =Tags(s)

MVALLEHYDHVPGDPEFTQLSHAVLLELVIDRLLLLLQSCMSYLENLGSEQMMPPAQAAGPCMSVGLTV
 R
 RFWDSLLRLGTLHQQLPQKGANQRETPTSKPTTKGEPARSPEYLTTKLIKPSSPVLGLPQTCQEPESIP
 VRASLQFPATTFKNTRSVHSQTIETALVPCDACASVQGSLQKVGVISLCQSQNLPSLGGFQQLVQDS
 MGLRPLPAATVGRWAAEQRKDLTRLSKHVEALRAQLEEAEGQKDGLRKQAGKLEQALKQEQGARRRQAE
 E
 DEQCLSEWEHDKQQLLTETSDLTKMATLERELKQQRESTQAVEAKAQQQLQEEGERRAAAERQVQQL
 Q
 VQQLEAQVQLLVGRLEGAGQQVCWASTELDKKARVDSMVRHQESLQAKQRALLKQLDQLDQEREELR
 GS
 LDEAEAQRARVEEQLQSEREQGQCQLRAQQELLQSLQREKQGLEQATTDLRLTILELERELEELKERERL
 LVAFPDLRPTETQIHGGRSSSVESQITCPTDSGNVTDHMERQVQANDIRIRVLQEENGRLQSMLSKIRE
 VAQQGGLKLIPQDRLWSPSSKGTQGATPPVQAKSTSPGPLGRQHLPSSRTGRTLLGQPCTSPPRQPCTSP
 PRQPCTSPPRQPCTSPSRQPCSQPSKSLLEGVTHLDTCTQNPIKVLVRLRKRLSPGRGQASSAHQPQERP
 M


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	78.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol


[View online »](#)

Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001017437</u>
Locus ID:	550631
UniProt ID:	<u>Q569K6</u>
RefSeq Size:	3023
Cytogenetics:	22q12.2
RefSeq ORF:	2103

Product images:



Coomassie blue staining of purified CCDC157 protein (Cat# [TP310351]). The protein was produced from HEK293T cells transfected with CCDC157 cDNA clone (Cat# [RC210351]) using MegaTran 2.0 (Cat# [TT210002]).