

Product datasheet for TP310363L

C1orf218 (NM_019049) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 1 open reading frame 218 (C1orf218), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210363 protein sequence Red=Cloning site Green=Tags(s)


MTKATPAVRTAYKFAKNHAKLGLKEVKSKLKHKENEEDYGTCSSTVQYTPVYKLHNEKGGNSEKRKLAQA
 RLKRPLKSLDGALYDDEDDDDIERASKLSSDGEASAYLYESDDSVETRVKTPYSGEMDLLGEILDTLS
 THSSDQGKLAAAKSLDFFRSMDDIDYKPTNKSNAPENNLAFLCGGSGDQAEWNLGQDDSAHGHKHL
 PS
 PRKRVSSSGLTDSLFIKEENS NKHLGADNVSDPTSLDFQLTSPEVSQTDKGKTEKRETLISQISDDLII
 PGLGRHSSTFVPWEKEGKEAKETSEDIGLLHEVVS LCHMTSDFQQSLNISDKNTNGNQ T

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	37.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_061922</u>
Locus ID:	54530


[View online »](#)

RefSeq Size: 2035
 Cytogenetics: 1q31.3
 RefSeq ORF: 1017
 Synonyms: FLJ20054, DKFZp547O0715

Product images:


Coomassie blue staining of purified C1orf218 protein (Cat# [TP310363]). The protein was produced from HEK293T cells transfected with C1orf218 cDNA clone (Cat# [RC210363]) using MegaTran 2.0 (Cat# [TT210002]).