

Product datasheet for TP310369


C4orf17 (NM_032149) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 4 open reading frame 17 (C4orf17), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210369 protein sequence Red =Cloning site Green =Tags(s) MNLNPPPTSALQIEGKGSIMARNVSCFLVRHTPHRRVCHIKGLNNIPICTVNDDENAFGTLWEVGQSNY LEKNRIPFANCYPPSTAVQKSPVRGMSPAPNGAKVPPRPHSEPSRKIKECFKTSSENPLVIKKEEIKAK RPPSPPKACSTPGSCSSGMTSTKNDVKANTICIPNYLDQEIKILAKLCSILHTDSLAEVLQWLLHATSKE KEWVSALIHSELAIEINLLTHHRRNTSMEPAAETGKPPTVKSPPTVKLPPNFTAKSKVLTRDTEGDQPTRV SSQGSEENKEVPKEAEHKPPLLIRNNMKIPVAEYFSKPNSPPRPNTQESGSAKPVSARSIQEYNLCPQR ACYPSTHRR TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	39.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_115525
Locus ID:	84103


[View online »](#)

UniProt ID: Q53FE4
 RefSeq Size: 1652
 Cytogenetics: 4q23
 RefSeq ORF: 1077

Product images:


Coomassie blue staining of purified C4orf17 protein (Cat# TP310369). The protein was produced from HEK293T cells transfected with C4orf17 cDNA clone (Cat# [RC210369]) using MegaTran 2.0 (Cat# [TT210002]).