

Product datasheet for TP310449

MRPL39 (NM_017446) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human mitochondrial ribosomal protein L39 (MRPL39), nuclear gene encoding mitochondrial protein, transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210449 protein sequence Red=Cloning site Green=Tags(s)

MEALAMGSRALRLWLVPAGGGIKWRFIATSPASQLSPTELTEMRNDFLNKEKARQLSLTPRTEKIEVKHV
GKTDPGTVFVMNKNISTPYSCAMHLSEWYCRKSILALVDGQPWDMYKPLTKSCEIKFLTFKDCDPGEVVK
AYWRSCAMMMGCVIERAFKDEYMNVLVRAPEVPVISGAFICYDVLDSKLDEWMPTKENLRSFTKDAHALI
YKDLPFETLEVEAKVALEIFQHSKYKVDFFIEEKASQNERIVKLRHIGDFIDVSEGPLIPRTSICFQYEV
SAVHNLQPTQPSLIRRFQGVSLPVHLRAHFTIWDKLLERSRKMVTEQSKATEECTST

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	38.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_059142
Locus ID:	54148



[View online »](#)

UniProt ID:	Q9NYK5
RefSeq Size:	1082
Cytogenetics:	21q21.3
RefSeq ORF:	1014
Synonyms:	C21orf92; L5mt; L39mt; MRP-L5; MRPL5; MSTP003; PRED22; PRED66; RPML5
Summary:	Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein. Two transcript variants encoding distinct isoforms have been described. A pseudogene corresponding to this gene is found on chromosome 5q. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified MRPL39 protein (Cat# TP310449). The protein was produced from HEK293T cells transfected with MRPL39 cDNA clone (Cat# [RC210449]) using MegaTran 2.0 (Cat# [TT210002]).