

## Product datasheet for **TP310458M**

### Nucleophosmin (NPM1) (NM\_001037738) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human nucleophosmin (nucleolar phosphoprotein B23, numatrin) (NPM1), transcript variant 3, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210458 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)  MEDSMDMDMSPLRPQNYLFGCELKADKDYHFKVDNDENEHQLSLRTVSLGAGAKDELHIVEAEAMNYEGS PIKVTLATLKMSVQPTVSLGGFEITPPVVLRLKCGSGPVHISGQHLVAVEEDAEESEEEEEVDKLLSISG KRSAPGGGSKVPQKKVKLADEDDDDDEDDDEDDDDDDDFDDEEAEEKAPVKKSIRDTPAKNAQKSNQN GKDSKPSSTPRSKGQESFKKQEKTPKTPKGPSSVEDIKAKMQASIEKAH  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	28.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_001032827</a></u>
Locus ID:	4869



[View online »](#)

UniProt ID: [P06748](#), [A0A140VJQ2](#)

RefSeq Size: 1347

Cytogenetics: 5q35.1


RefSeq ORF: 777

Synonyms: B23; NPM

**Summary:** The protein encoded by this gene is involved in several cellular processes, including centrosome duplication, protein chaperoning, and cell proliferation. The encoded phosphoprotein shuttles between the nucleolus, nucleus, and cytoplasm, chaperoning ribosomal proteins and core histones from the nucleus to the cytoplasm. This protein is also known to sequester the tumor suppressor ARF in the nucleolus, protecting it from degradation until it is needed. Mutations in this gene are associated with acute myeloid leukemia. Dozens of pseudogenes of this gene have been identified. [provided by RefSeq, Aug 2017]

**Protein Families:** Druggable Genome, Stem cell - Pluripotency, Transcription Factors

### Product images:



Coomassie blue staining of purified NPM1 protein (Cat# [TP310458]). The protein was produced from HEK293T cells transfected with NPM1 cDNA clone (Cat# [RC210458]) using MegaTran 2.0 (Cat# [TT210002]).