

Product datasheet for TP310514M

Plakophilin 2 (PKP2) (NM_001005242) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of Human plakophilin 2 (PKP2), transcript variant 2a, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210514 protein sequence Red =Cloning site Green =Tags(s)

MAAPGAPAEYGYIRTVLGQQILGQLDSSSLALPSEAKLKLKAGSSGRGGQTVKSLRIQEQVQQLARKGRS
SVGNLHRTSSVPEYVYNLHLVENDFVGGRRSPVPKTYDMLKAGTTATYEGRWGRGTAQYSSQKSVEERS
LRHPLRRLEISPDSSPERAHYTHSDYQYSQRSQAGHTLHHQESRRAALLVPPRYARSEIVGVSRAQTTSR
QRHFDYHRQYQHGSVSDTVFDSIPANPALLTYPRPGTSRSMGNLLEKENYLTAGLTVGQVRPLVPLQPV
TQNRASRSSWHQSSFHSTRTLREAGPSVAVDSSGRRRAHLTVGQAAAGGSGNLLTERSTFTDSQLGNADME
MTLERAVSMLEADHMPPSRISAAATFIQHECFQKSEARKRVNQLRGILKLLQLLKVQNEVDQRAVCGALR
NLVFEDNDNKLEVAELNGVPRLLQVLKQTRDLETKKQITGLLWNLSSNDKLNLMITEALLTLTENIIIP
FSGWPEGDYPKANGLLDFIFYNVTGCLRNMSSAGADGRKAMRRCDGLIDSLVHYVRGTIADYQPPDDKAT
ENCVCILHNLSYQLEAELPEKYSQNIYIQRNRIQTDNNSIGCFGSRSRKVKQEYQDVPMPPEEKSNPKGV
EWLWHSIVIRMYLSLIAKSVRNYTQEASLGALQNLTAGSGPMPTSAQTVMQKESGLQHTRKMLHVGDPD
VKKTAISLLRNLRLSLQNEIAKETLPDLVSIIPDTPSTDLLIETTASACYTLNIIQNSYQNARDLL
NTGGIQKIMASAGDAYASNKASKAASVLLYSLWAHTELHHAYKKAQFKKTDVNSRTAKAYHSLKD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	92.7 kDa
Concentration:	>0.05 µg/µL as determined by Bradford protein assay method.
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	100 mM Glycine, pH 3.5, 10% Glycerol
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_064303.2




[View online »](#)

Locus ID: 5318
UniProt ID: [Q99959](#)
RefSeq Size: 4307
Cytogenetics: 12p11.21
RefSeq ORF: 2511
Synonyms: ARVD9

Summary: This gene encodes a member of the arm-repeat (armadillo) and plakophilin gene families. Plakophilin proteins contain numerous armadillo repeats, localize to cell desmosomes and nuclei, and participate in linking cadherins to intermediate filaments in the cytoskeleton. This gene product may regulate the signaling activity of beta-catenin. Two alternately spliced transcripts encoding two protein isoforms have been identified. A processed pseudogene with high similarity to this locus has been mapped to chromosome 12p13. [provided by RefSeq, Jul 2008]

Protein Pathways: Arrhythmogenic right ventricular cardiomyopathy (ARVC)

Product images:



Coomassie blue staining of purified PKP2 protein (Cat #[TP310514]). The protein was produced from mammalian cells transfected with PKP2 cDNA clone (Cat #[RC210514]).