

Product datasheet for TP310515L

GSDME (NM_004403) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human deafness, autosomal dominant 5 (DFNA5), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210515 protein sequence Red =Cloning site Green =Tags(s)

MFAKATRNFLREVDADGDLIAVSNLNDSKQLLSLVTKKKRFWCWQRPKYQFLSLTLGDVLIEDQFPSP
 VVVESDFVKYEGKFANHVSGTLETALGKVKLNLLGGSSRVESQSSFGTLRKQEVLDLQQLIRDSAERTINLR
 NPVLQQVLEGRNEVLCVLTQKITTMQKCVISEHMQVEEKCGGIVGIQTKTVQVSATEDGNVTKDSNVVLE
 IPAATTIAYGVIELYVKLDGLFFECLLRGKQGGFENKKRIDSVYLDHLVREFAFIDMPDAAHGISSQDG
 PLSVLKQATLLLRNFHPFAELPEPQQTALSDIFQAVVFDDELLMVLEPVCDDLVSGLSPTVAVLGELKP
 RQQQDLVAFLQLVGCSLQGGCPGEDAGSKQLFMTAYFLVSALAEMPDSAAALLGTCCCKLQIIPTLCHLL
 RALSDDGVSDLEDPTLTPLKDTERFIVQRLFASADISLERLKS SVKAVILKDSKVFLLLCITLNLGLCA
 LGREHS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	54.4 kDa
Concentration:	>0.1 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_004394
Locus ID:	1687
UniProt ID:	O60443 , A0A024RA58
RefSeq Size:	2521
Cytogenetics:	7p15.3
RefSeq ORF:	1488
Synonyms:	DFNA5; ICERE-1
Summary:	Hearing impairment is a heterogeneous condition with over 40 loci described. The protein encoded by this gene is expressed in fetal cochlea, however, its function is not known. Nonsyndromic hearing impairment is associated with a mutation in this gene. Three transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome

Product images:

Coomassie blue staining of purified GSDME protein (Cat# [TP310515]). The protein was produced from HEK293T cells transfected with GSDME cDNA clone (Cat# [RC210515]) using MegaTran 2.0 (Cat# [TT210002]).