

Product datasheet for **TP310524L**

UBXN2A (NM_181713) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human UBX domain protein 2A (UBXN2A), 1 mg

Species: Human

Expression Host: HEK293T

**Expression cDNA Clone
or AA Sequence:** >RC210524 protein sequence
Red=Cloning site **Green**=Tags(s)

MKDVDNLKSIKEEWVCETGSDNQPLGNNQQSNCEYFVDSLFEAAQKVSSKCVSPAEEQKKQVDVNIKWLWKN
GFTVNDDFRSYSDGASQQFLNSIKKGELPSELQGIFDKEEVDVKVEDKKNEICLSTKPVFQPFSGQGHRL
GSATPKIVSKAKNIEVENKNNLSAVPLNNLEPITNIQIWLANKRIVQKFNITHRVSHIKDFIEKYQGSQ
RSPPFSLATALPVLRLLEDLTLLEADLQNAVIIQRLQKTASFRELSEH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 29.1 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_859064](#)


Locus ID: 165324

UniProt ID: [P68543](#)



[View online »](#)

RefSeq Size: 1909
Cytogenetics: 2p23.3
RefSeq ORF: 777
Synonyms: UBXD4

Product images:

Coomassie blue staining of purified UBXN2A protein (Cat# [TP310524]). The protein was produced from HEK293T cells transfected with UBXN2A cDNA clone (Cat# [RC210524]) using MegaTran 2.0 (Cat# [TT210002]).