

Product datasheet for **TP310533L**

ER81 (ETV1) (NM_004956) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ets variant 1 (ETV1), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210533 protein sequence Red =Cloning site Green =Tags(s)
	<p>MDGFYDQQVPYMTNSQRGRNCNEKPTNVRKRKFINRDLAHDSEELFQDLSQLQETWLAEAQVPDNDEQF VPDYQAESLAFHGLPLKIKKEPHSPCSEISSACSQEQPFKFSYGEKCLYNVSAYDQKPQVGMRPSNPPTP SSTPVSPHLLHASPNSHTPKPDRAFPAPHLPPSQSIPDSSYPMDHRFRRLSEPCNSFPPLPTMPREGRPM YQRQMSEPNIFFPPQGFKQEYHDPVYEHNTMVGSAASQSFPPPLMIKQEPDRDFAYDSEVPSCHSIYMRQE GFLAHPSTRTEGCMFEKGRQFYDDTCVPEKFDGDIKQEPGMYREGPTYQRRGSLQLWQFLVALLDDPSN SHFIAWTGRGMEFKLIEPEEVARRWGIQKNRPAMNYDKLSRSLRYYYEKGIMQKVAGERYVYKFCVCDPEA LFSMAFPDNQRPLLKTDMERHINEEDTVPLSHFDESMAYMPEGGCCNPHYPYNEGYVY</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	55 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_004947</u>




[View online »](#)

Locus ID: 2115
UniProt ID: [P50549](#)
RefSeq Size: 6824
Cytogenetics: 7p21.2
RefSeq ORF: 1431
Synonyms: ER81

Summary: This gene encodes a member of the ETS (E twenty-six) family of transcription factors. The ETS proteins regulate many target genes that modulate biological processes like cell growth, angiogenesis, migration, proliferation and differentiation. All ETS proteins contain an ETS DNA-binding domain that binds to DNA sequences containing the consensus 5'-CGGA[AT]-3'. The protein encoded by this gene contains a conserved short acidic transactivation domain (TAD) in the N-terminal region, in addition to the ETS DNA-binding domain in the C-terminal region. This gene is involved in chromosomal translocations, which result in multiple fusion proteins including EWS-ETV1 in Ewing sarcoma and at least 10 ETV1 partners (see PMID: 19657377, Table 1) in prostate cancer. In addition to chromosomal rearrangement, this gene is overexpressed in prostate cancer, melanoma and gastrointestinal stromal tumor. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2016]

Protein Families: ES Cell Differentiation/IPS, Transcription Factors

Product images:



Coomassie blue staining of purified ETV1 protein (Cat# [TP310533]). The protein was produced from HEK293T cells transfected with ETV1 cDNA clone (Cat# [RC210533]) using MegaTran 2.0 (Cat# [TT210002]).